

AD-A227 924

THE DEVELOPMENT OF A SYSTEM FOR
ENHANCED RADIONUCLIDE DETECTION

ANNUAL REPORT

JOHN K. FROST
KAREN E. WARDEN
ISAAC N. BANKMAN
NORMAN J. PRESSMAN
PRABODH K. GUPTA
DARRELL L. COOK
GARY W. GILL
RICHARD L. SHOWERS
RICHARD K. TRAUB

NOVEMBER 1987

DTIC
ELECTE
OCT 23 1990
S B D
Ca

Supported by

U.S. ARMY MEDICAL RESEARCH AND DEVELOPMENT COMMAND
Fort Detrick, Frederick, Maryland 21701-5012

Contract No. DAMD17-84-C-4222

The Johns Hopkins University School of Medicine
Baltimore, Maryland 21205

Approved for public release; distribution unlimited.

The findings in this report are not to be construed as an
official Department of the Army Position unless so
designated by other authorized documents.


90 10 22 1989

REPORT DOCUMENTATION PAGE				Form Approved OMB No. 0704-0188	
1a. REPORT SECURITY CLASSIFICATION Unclassified			1b. RESTRICTIVE MARKINGS		
2a. SECURITY CLASSIFICATION AUTHORITY			3. DISTRIBUTION/AVAILABILITY OF REPORT Approved for public release; distribution unlimited.		
2b. DECLASSIFICATION/DOWNGRADING SCHEDULE					
4. PERFORMING ORGANIZATION REPORT NUMBER(S)			5. MONITORING ORGANIZATION REPORT NUMBER(S)		
6a. NAME OF PERFORMING ORGANIZATION The Johns Hopkins University School of Medicine		6b. OFFICE SYMBOL (if applicable)	7a. NAME OF MONITORING ORGANIZATION		
6c. ADDRESS (City, State, and ZIP Code) Baltimore, Maryland 21205			7b. ADDRESS (City, State, and ZIP Code)		
8a. NAME OF FUNDING/SPONSORING ORGANIZATION U.S. Army Medical Research & Development Command		8b. OFFICE SYMBOL (if applicable)	9. PROCUREMENT INSTRUMENT IDENTIFICATION NUMBER Contract No. DAMD17-84-C-4222		
8c. ADDRESS (City, State, and ZIP Code) Fort Detrick Frederick, Maryland 21701-5012			10. SOURCE OF FUNDING NUMBERS		
			PROGRAM ELEMENT NO. 63764A	PROJECT NO. 3M4- 63764D995	TASK NO. AB
11. TITLE (Include Security Classification) (U) The Development of a System for Enhanced Radionuclide Detection					
12. PERSONAL AUTHOR(S) John K. Frost, Karen E. Warden, John K. Frost, Jr., Isaac N. Bankman, Norman J. Pressman, Prabodh K. Gupta, Darrell L. Cook, G.W. Gill, R.L. Showers, R.K. Traub					
13a. TYPE OF REPORT Annual		13b. TIME COVERED FROM 8/1/85 TO 7/31/86		14. DATE OF REPORT (Year, Month, Day) 1987 November	
15. PAGE COUNT					
16. SUPPLEMENTARY NOTATION					
17. COSATI CODES			18. SUBJECT TERMS (Continue on reverse if necessary and identify by block number) Autoradiography; Computer; Detection; Energy Transfer, Fluorochrome; Histology; Image Analysis; Image Intensifica- tion; Instrumentation; Photonic, Photo-limited Imaging		
FIELD	GROUP	SUB-GROUP			
19. ABSTRACT (Continue on reverse if necessary and identify by block number) Frequently, toxicologic processes involve concentrations of biologically active substances so low that to be detected, localized, and quantified they require attachment of appropriate radiolabels and subsequent autoradiography exposure periods that may extend for many months. The major objective of this research is to develop and deliver a system to reduce this long data collection period while maintaining the advantage of spatial resolution inherent in autoradiography. After completing concept validation and after careful testing of three low-light-level camera system which were submitted by vendors for evaluation, a photon-limited camera and imaging system by Hamamatsu and a programmable image display unit by Gould were selected and ordered for purchase. Software programs to distinguish and analyze microstructures and their properties were generated.					
20. DISTRIBUTION/AVAILABILITY OF ABSTRACT <input type="checkbox"/> UNCLASSIFIED/UNLIMITED <input checked="" type="checkbox"/> SAME AS RPT <input type="checkbox"/> DTIC USERS			21. ABSTRACT SECURITY CLASSIFICATION Unclassified		
22a. NAME OF RESPONSIBLE INDIVIDUAL Mary Frances Bostian			22b. TELEPHONE (Include Area Code) 301-663-7325		22c. OFFICE SYMBOL SGRD-RMI-S

SUMMARY

Frequently, toxicologic processes involve concentrations of biologically active substances so low that to be detected, localized, and quantified they require attachment of appropriate radiolabels and subsequent autoradiography exposure periods that may extend for many months. The major objective of this research is to develop and deliver a system to reduce this long data collection period while maintaining the advantage of spatial resolution inherent in autoradiography.

After completing concept validation and after careful testing of three low light level camera systems which were submitted by vendors for evaluation, a photon-limited camera and imaging system by Hamamatsu and a programmable image display unit by Gould were selected and ordered for purchase. Software programs to distinguish and analyze microstructures and their properties were generated.



Accession For	
NTIS GRA&I	<input checked="" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By	
Distribution/	
Availability Codes	
Dist	Avail and/or Special
A-1	

FOREWORD

Citations of commercial organizations and trade names in this report do not constitute an official Department of the Army endorsement or approval of the products or services of these organizations.

TABLE OF CONTENTS

SUMMARY	1
FOREWORD	2
I. INTRODUCTION	5
II. MATERIALS AND METHODS	6
A. Data Acquisition and Processing Systems	6
B. Development of Computer Software for Data Analysis	6
III. RESULTS	8
A. Data Acquisition and Processing Systems	8
1. Gould IP8400 Programmable Image Display Unit	8
2. Hamamatsu C1966-20 Video-Intensified Microscopy System (VIM)	8
B. Development of Computer Software for Data Analysis	10
C. Literature Review	19
IV. DISCUSSION	20
V. CONCLUSIONS	22
BIBLIOGRAPHY	23
GLOSSARY	31
APPENDICES	33
A. MSAR Flowcharts	33
B. MSAR Program Code	43

LIST OF FIGURES

Figure 1.	Concept Validation Image #ED0030 for Signal Dosimetry Experiment with One-second Summation Time	13
Figure 2.	Microstructural Analysis Routine Test Image	16
Figure 3.	Test Patterns of Pixel Contiguousness for Different Threshold Ranges of the Microstructure Detection Test Image	17
Figure 4.	Plot of Image Area above Threshold versus Signal Integration Time	21

LIST OF TABLES

Table 1.	Tabulated Results from the Particlet Processing of Experimental Image #ED0030	14
Table 2.	Tabulated Results from the MSAR Processing of Experimental Image #ED0030	14
Table 3.	Tabulated Results for a Particular Threshold Range from MSAR Processing of Experimental Image #ED0030 (Abridged)	15
Table 4.	Tabulated Results from MSAR Processing of Microstructure Detection Test Image	18

I. INTRODUCTION

The U.S. Army Medical Research and Development Command has a need to measure and accurately localize attogram concentrations of chemical compounds in animal and human tissues. Indirect enhancement of the specific energies of the molecules of these compounds is a fundamental approach to this problem.⁵⁹

Enhancement involves the substitution of elements in the compound with their radioactive isotopes during synthesis. Less energetic radioactive isotopes, such as ^{14}C (0.156 MeV) and ^3H (0.01861 MeV), do not present difficult specimen handling problems and have half-lives of 5,570 and 12.3 years, respectively. They permit well-designed biological experiments to be conducted. The observed mode of decay for both ^{14}C and ^3H is β -emission. Depending upon the isotope being used, current methods require contact autoradiographic exposures of as much as 6 months to one year duration to record a usable image.⁹⁰

Scintillation detection and counting by this method will be a more effective alternative, for detecting weak β -emitting radionuclides, than autoradiography. In this method, the energy of the β -particles excites a scintillator, which converts β -particle kinetic energy into potentially detectable photons. In this current effort, these photonic events are temporally and spatially accumulated as a function of position. The problem is being approached by demonstrating the feasibility of rapidly and microscopically imaging low-concentration radiolabeled tissues by developing: (1) methods for converting the energy of emitted β -particles to photons and (2) working instrumentation for recording the resulting scintillations with high spatial and temporal resolution.^{12, 23, 31, 60, 91, 97}

II. MATERIALS AND METHODS

A. Data Acquisition and Processing Systems

There are two major items of equipment whose purchases are essential for this contractual effort:

- 1) The Gould IP8400 Programmable Image Display Unit*, and
- 2) The Hamamatsu C1966-20 Video-Intensified Microscopy System (VIM)**.

While much of the thorough evaluation of each of these items was carried out during the previous year, it was completed in this year with decision to purchase. A request for approval was granted, and the order for this equipment was placed May 23, 1986.

B. Development of Computer Software for Data Analysis

The program Particlet was obtained from our USAMRICD collaborator, R.K. Traub. The primary function of Particlet is to enumerate and to calculate the area of microstructures having intensity values above or below a user-supplied threshold. The results displayed to the user are:

- 1) threshold value,
- 2) particle count, and
- 3) total integrated area of counted particles.

In this context, a "particle" is defined as a contiguous set of pixels whose values all exceed or equal some threshold, and which are surrounded entirely by pixels whose values are all below that threshold.

Two modifications had to be made before Particlet could be incorporated onto our DEC VAX 11/780 system, worked upon to develop the data analysis program, and delivered to the Army in the format required on completion of the contract:

- 1) The first required the alteration of program code, so that image data were read in the form of byte-packed 8-bit data and then converted to the program format of 64-bit data.
- 2) The second was to increase the system working disk space to fifty thousand pages, because of the large virtual memory.

* Gould, Inc., San Jose, CA

** Photonic Microscopy, Inc., Oakbrook, IL

A program called Microstructural Analysis Routine (MSAR) was then developed for this project by our computer scientist. This new program uses Particlet as a basic concept model and is written to analyze data for which Particlet had not been designed. MSAR thus offers a different data analysis and output of the extremely limited photon events (Section III. B.).

III. RESULTS

A. Data Acquisition and Processing Systems

1. Gould IP8400 Programmable Image Display Unit

The Gould IP8400 is a programmable digital image array processing and display unit. It contains look-up tables and special hardware processors to allow user manipulation of stored images and Random Access Memory (RAM) for image storage. In its present configuration it is attached to a VAX-11/780, as its host computer, as a peripheral device on the DEC Unibus.

Image data, acquired from the Hamamatsu C1966 VIM processor by the VAX-11/780, are transferred to the IP8400 RAM at a minimum rate of 1.5 Mbytes/sec. Currently the IP8400 has over 524 Kbytes of RAM in two 512 x 512 x 8-bit memory channels and in look-up tables and registers. Two users can access the IP8400 concurrently, but memory cannot be shared. Each display channel consists of one memory channel interfaced to a Video Output Controller (VOC) capable of driving a monochrome monitor or an RGB color monitor with composite video signals.

The master timing of the IP8400 is selectable and can be changed by software calls to internal registers. Each display channel can be manipulated by an associated joystick and has an alphanumeric overlay generator for producing annotation on the display monitor.

The IP8400 is also equipped with a high-speed pipeline processor, the Digital Video Processor (DVP), which will operate on a combination of 512 x 512 x 8-bit memories. It allows real-time operations such as multiply, add, subtract, compare, and reciprocal divide.

2. Hamamatsu C1966-20 Video-Intensified Microscopy System (VIM)

This is an extremely sensitive system developed for use in ultra low-light-level situations. Its characteristics include a wide dynamic range (to image over different intensity levels), high spatial resolution, a high signal-to-noise ratio, and real time data acquisition capabilities.

After evaluating three low-light-level cameras in the first year (Zeiss Venus, Dage-MTI, and the Hamamatsu C1966-20) and final system evaluations in this second year, the Hamamatsu C1966-20 VIM was chosen for purchase on the basis of its sensitivity and photon-limited characteristics.

The system converts emitted photons into a proportional number of electrons through the use of a two-dimensional photon counting tube that consists of three parts. On one face is a photocathode, where the initial optical image is focused. The photons are converted to electrons here and

subsequently accelerated to the second part of the counting tube, a two-stage microchannel plate (MCP), where they are amplified by electron avalanching. The MCP has two important qualities:

- a) due to its design, the spatial integrity of the image is maintained and
- b) because there is an adjustable gain setting, image intensifier gain (avalanching) can be controlled by the user.

The gain of the system can be operated in an analog mode, if sufficient light is present, or in digital mode for low-light-level imaging. Once the amplified electrons strike the phosphor screen that follows the MCP, a resultant image is presented to an optically coupled low lag vidicon. This produces signal pulses which may then be discriminated, with an acceptance threshold, to eliminate thermal noise that may occur.

The resolution approaches that of a Silicon Intensified Target (SIT) camera (300 to 350 TV lines as compared with 400 for the SIT). By selecting the type of intensifier and associated photocathode, a range of spectral sensitivity from 140 to 850 nm is possible.

The following is a list of the selected components:

- a) Specially selected ultra-low noise/high resolution VIM photon counting camera head with optional Peltier cooling system. The Peltier cooling system was evaluated to determine its effectiveness in minimization of thermal noise. In addition, vendor specifications indicate that a dark count range of 5 to 10 cps can be achieved by carefully selecting the image intensifier tube.
- b) VIM Controller and control keyboard. To reduce the number of incident light sources and thereby assure a light-tight laboratory environment, the keyboard has been adapted with a disable switch.
- c) DEC DMA (direct memory access) interface and software. This feature directly couples the C1966-20 to the DEC VAX 11/780 computer memory, allowing digitized image data or control/command language to be transferred in either direction between the two devices. It offers the advantage of reducing the amount of user interaction necessary for operation of routine functions.

B. Development of Computer Software for Data Analysis

The Microstructural Analysis Routine (MSAR) was developed to enhance the performance of the USAMRICD Program Particlelet (Section II.B.), by providing additional functionality, in order to meet the need to quantitatively assess images derived from radionuclide-labeled sections obtained by the System under development. MSAR (Appendix A and Appendix B) has the ability to perform the following five functions:

1. Read Image Data into VAX 11/780 Memory
2. Define Region of Interest (ROI)
3. Select Threshold
4. Determine Adjacency Criteria
5. Calculate Microstructure Statistics

1. Read Image Data into VAX 11/780 Memory

The user may read a monochromatic, 512 x 512 pixel image from disk into VAX memory. The user is also requested to supply the physical dimension of the pixel in microns, if it is known. This factor is important in calculating the areas and perimeters of the microstructures.

2. Define Region of Interest (ROI)

A region of interest (ROI), which is defined by the user, is a subset of the full field image on which microstructure analysis is performed. If no region of interest is specified, the entire image field is used. The user selects an ROI by entering the X, Y coordinate values for the upper left and lower right corners of the region.

3. Select Threshold

Three types of thresholding specifications are available for user selection:

- a) Minimum Threshold — all pixel values BELOW this value will be ignored during analysis,
- b) Maximum Threshold — all pixel values ABOVE this value will be ignored during analysis,
- c) Threshold Range — a combination of the above specifications such that only pixel values lying BETWEEN the minimum AND maximum, inclusively, will be used in analysis.

These selected thresholds must be values between 0 and 255 to be valid. In addition, the user is given the option of specifying a threshold increment for both the minimum and maximum threshold, either individually or together. This allows analysis to be performed in sequential steps automatically, thereby eliminating user interaction for repetitive operations.

4. Determine Adjacency Criteria

The user is given two choices in determining the criteria for contiguity of adjacent points: 1. Orthogonal; 2. Orthogonal and diagonal. If no adjacency mode is selected, the default is to number two.

5. Calculate Microstructure Statistics

MSAR performs a variety of microstructure calculations. For each specified threshold range of an image, the following are calculated:

- a) Total area of all microstructures (in pixels)
- b) Number of microstructures
- c) Average object area
- d) Average X,Y center of region of interest
- e) Average X,Y weighted center of region of interest
- f) Total number of perimeter points
- g) Average number of perimeter points per object
- h) Total perimeter length
- i) Average perimeter length
- j) Average horizontal and vertical feret diameters
- k) Cumulative and average (per object) integrated optical densities
- l) Average maximum optical density (per object)

The perimeter is calculated in two ways. The first, the total number of perimeter points, is the sum of all of the points on the boundary of the microstructure. The second, the total perimeter length, is the distance travelled, in pixels (or microns if the user supplied a pixel to micron conversion factor), of one complete trace of the perimeter. Once the microstructures have been uniquely labeled, the following statistical calculations are performed for each microstructure within a particular threshold range:

- a) X,Y location of the geometric center
- b) X,Y weighted center
- c) Area
- d) Total number of perimeter points and length
- e) Maximum horizontal and vertical diameters
- f) Integrated optical density (pixel values)
- g) Average integrated optical density (pixel values)
- h) Maximum optical density (pixel values)

Many of the calculation results from the preceding lists are given in units of length. For this reason, the pixel size is requested at the onset of the program. The results obtained will be more meaningful when MSAR is allowed to make calculations using actual dimensions.

The hardcopy results obtained upon the completion of Particlet and MSAR processing of the Experimental Image #ED0030 (Figure 1.) are given in Tables 1 and 2. While each routine offers the number of microstructures and the average area of these within a threshold range (or "slice level"), it is evident from this data that the results from MSAR are more extensive and useful. The information concerning averages is helpful in determining general data trends. For instance, the average weighted X and Y center can direct the user to the region of highest intensity, which corresponds to the region with the highest concentration of radiolabeled substance.

An additional level of detail in MSAR is illustrated in Table 3. Here information about specific microstructures is available for each unique microstructure in the image.

FIGURE 1. Concept Validation Image #ED0030 for Signal Dosimetry
Experiment with One-second Summation Time.




TABLE 1. Tabulated Results from the Particlet Processing of
Experimental Image #ED0030

ED0030. DAT

slice level,	number,	area of particles
10.00000	1	248303
20.00000	345	1246
30.00000	16	47
40.00000	4	15
50.00000	2	12
60.00000	2	12
70.00000	2	12
80.00000	2	11
90.00000	2	10
100.0000	2	10
110.0000	2	8
120.0000	2	6
130.0000	2	6
140.0000	2	5
150.0000	2	5
160.0000	2	4
170.0000	2	4
180.0000	2	4
190.0000	2	2
200.0000	1	1
210.0000	0	0
220.0000	0	0
230.0000	0	0
240.0000	0	0
250.0000	0	0

TABLE 2. Tabulated Results from the MSAR Processing of
Experimental Image #ED0030

FILENAME = ED0030.DAT

REGION OF INTEREST - TOP LEFT (X= 1, Y= 1), BOTTOM RIGHT (X=512, Y=512)

THRESHOLD RANGE	TOTAL AREA	NUMBER OBJECTS (MICROSTRUCTURES)	AVERAGE OBJECT AREA	AVERAGE X- CENTER	AVERAGE Y- CENTER	AVERAGE WEIGHTED X- CENTER	AVERAGE WEIGHTED Y- CENTER
20-255	1246	345	3.61	255.19	303.44	234.48	309.66

THRESHOLD RANGE	TOTAL PERIMETER POINTS	AVERAGE PERIMETER POINTS PER OBJECT	TOTAL PERIMETER	AVERAGE PERIMETER PER OBJECT	AVERAGE HORIZ FERET DIAMETER	AVERAGE VERTICAL FERET DIAMETER	CUMULATIVE 100	AVERAGE 100 PER OBJECT	AVERAGE ODD=1 PER OBJECT
20-255	287	0.83	346.00	1.00	2.44	1.45	31036	89.96	25.80

TABLE 3. Tabulated Results for a Particular Threshold Range from
MSAR Processing of Experimental Image #ED0030 (Abridged)

FILENAME = ED0030_020_255.MBT LOWER THRESHOLD = 20 UPPER THRESHOLD = 255
 345 UNIQUE MICROSTRUCTURES FOUND
 DEFAULT PIXEL SIZE = 1.000 MICRONS FLAG = 1

UNIQUE MICROSTRUCTURE	X CENTER	Y CENTER	WEIGHTED X CENTER	WEIGHTED Y CENTER	AREA	PERIMETER COUNT	PERIMETER LENGTH	MAXIMUM HORIZONTAL RUN-LENGTH	MAXIMUM VERTICAL RUN-LENGTH	AVERAGE 100	AVERAGE 100	QDmax
1	223.00	27.43	223.94	27.43	8.00	-1	-1.00	3	2	212	26.50	30
2	223.00	30.00	222.94	30.00	3.00	-1	-1.00	3	1	75	23.00	28
3	54.30	37.00	54.49	37.00	2.00	-1	-1.00	2	1	51	24.50	24
4	54.00	39.00	53.49	39.00	3.00	-1	-1.00	3	1	77	22.47	28
5	108.00	39.00	108.00	39.00	1.00	-1	-1.00	1	1	23	23.00	23
6	220.30	39.00	220.47	39.00	2.00	-1	-1.00	2	1	43	21.50	23
7	323.30	42.00	323.56	42.00	2.00	-1	-1.00	2	1	48	24.00	27
8	108.00	43.00	107.90	43.00	3.00	-1	-1.00	3	1	67	29.00	33
9	243.00	43.00	242.95	43.00	2.00	-1	-1.00	2	1	77	22.47	27
10	174.30	44.00	174.36	44.00	2.00	-1	-1.00	2	1	48	24.00	27
11	164.00	43.00	163.99	43.00	3.00	-1	-1.00	3	1	74	21.33	27
12	228.30	43.00	228.30	43.00	2.00	-1	-1.00	2	1	44	23.00	23
13	138.30	46.00	138.44	46.00	2.00	-1	-1.00	2	1	43	22.50	23
14	206.30	49.00	206.49	49.00	4.00	-1	-1.00	4	1	97	24.25	26
15	248.30	72.00	248.48	72.00	2.00	-1	-1.00	2	1	44	23.00	24
16	249.00	74.00	248.97	74.00	3.00	-1	-1.00	3	1	72	24.00	26
17	434.30	74.00	434.43	74.00	2.00	-1	-1.00	2	1	44	23.00	24
18	144.30	79.00	144.33	79.00	2.00	-1	-1.00	2	1	43	22.50	24
19	290.00	79.00	289.94	79.00	3.00	-1	-1.00	3	1	70	23.33	23
20	237.30	80.00	237.44	80.00	2.00	-1	-1.00	2	1	43	22.50	23
21	290.00	81.00	289.94	81.00	3.00	-1	-1.00	3	1	70	23.33	23
22	348.30	82.00	348.47	82.00	10.00	-1	-1.00	4	3	240	24.00	26
23	212.30	89.00	212.33	89.00	2.00	-1	-1.00	2	1	45	22.50	24
24	229.00	92.00	228.93	92.00	3.00	-1	-1.00	3	1	70	23.33	23
25	182.00	97.00	182.00	97.00	1.00	-1	-1.00	1	1	20	20.00	20
26	344.00	99.00	344.00	99.00	1.00	-1	-1.00	1	1	21	21.00	21
27	336.00	103.00	335.93	103.00	3.00	-1	-1.00	3	1	70	23.33	23
28	326.30	106.00	326.49	106.00	2.00	-1	-1.00	2	1	47	23.50	24
29	323.30	110.00	323.36	110.00	2.00	-1	-1.00	2	1	45	22.50	23
30	417.00	111.00	417.00	111.00	1.00	-1	-1.00	1	1	20	20.00	20
31	243.30	112.00	243.46	112.00	2.00	-1	-1.00	2	1	44	22.00	22
32	417.30	113.00	417.24	113.00	2.00	-1	-1.00	2	1	44	23.00	23
33	246.00	114.00	245.96	114.00	3.00	-1	-1.00	3	1	49	23.00	23
34	290.00	121.00	289.94	121.00	3.00	-1	-1.00	3	1	49	23.00	23
35	280.30	123.00	280.50	123.00	2.00	-1	-1.00	2	1	50	21.00	23
36	29.30	124.00	29.49	124.00	2.00	-1	-1.00	2	1	49	24.50	23
37	89.00	127.00	89.02	127.00	2.00	-1	-1.00	2	1	81	27.00	31
38	28.30	128.00	28.36	128.00	2.00	-1	-1.00	2	1	45	22.50	23
39	88.30	129.00	88.46	129.00	4.00	-1	-1.00	4	1	105	33.75	38
40	336.40	133.10	336.38	133.04	10.00	-1	-1.00	4	4	211	23.40	26
41	91.17	137.08	91.10	137.17	12.00	11	14.00	5	3	312	26.00	37
42	358.30	142.00	358.31	142.00	2.00	-1	-1.00	2	1	47	23.50	24
43	387.00	142.00	387.00	142.00	1.00	-1	-1.00	1	1	24	26.00	26
44	90.30	144.00	90.45	144.00	2.00	-1	-1.00	2	1	44	22.00	24
45	89.00	146.00	89.00	146.00	1.00	-1	-1.00	1	1	25	22.00	23


To examine the integrity and accuracy of the Microstructural Analysis Routine, a test image was synthesized (Figure 2).

FIGURE 2. Microstructural Analysis Routine Test Image

Figure 2a. Digital Representation

10	10	20	20	0	30
0	20	50	40	10	0
10	20	10	30	10	20
0	10	10	0	40	0
0	30	30	40	50	30
30	20	20	0	40	40

FIGURE 2b. Analog Representation from Monitor



This test image consists of a six by six matrix containing multiples of six different numerical (pixel) values ranging from zero to fifty in increments of ten. These values were arranged so that, as MSAR processed the image (first by analyzing the 0 to 255 threshold range, then 10 to 255, and so on, until the lower threshold value was equal to fifty), the pixels that fell within the threshold range formed test patterns representing the various conditions of pixel contiguousness (Figure 3). Table 4 shows the results of MSAR execution on the test matrix image.

FIGURE 3. Test Patterns of Pixel Contiguosness for Different Threshold Ranges of the Microstructure Detection Test Image




TABLE 4. Tabulated Results from MSAR Processing of
Microstructure Detection Test Image

Table 4a

FILENAME = TSTOBJ.DAT

REGION OF INTEREST = TOP LEFT (X= 1, Y= 1), BOTTOM RIGHT (X= 250, Y= 250)

THRESHOLD RANGE	TOTAL AREA	NUMBER OF OBJECTS (MICROSTRUCTURES)	AVERAGE OBJECT AREA	AVERAGE X-CENTER	AVERAGE Y-CENTER	AVERAGE WEIGHTED X-CENTER	AVERAGE WEIGHTED Y-CENTER
0-255	26	1	26.00	3.30	3.30	3.30	3.30
10-255	28	2	14.00	4.70	3.30	3.30	3.30
20-255	20	4	5.00	4.74	2.84	2.84	2.84
30-255	12	4	3.25	3.79	3.41	3.41	3.41
40-255	7	2	3.50	4.23	3.40	3.40	3.40
50-255	2	2	1.00	4.00	3.50	3.50	3.50
60-255	0	0	0.00	0.00	0.00	0.00	0.00
70-255	0	0	0.00	0.00	0.00	0.00	0.00
80-255	0	0	0.00	0.00	0.00	0.00	0.00
90-255	0	0	0.00	0.00	0.00	0.00	0.00
100-255	0	0	0.00	0.00	0.00	0.00	0.00
110-255	0	0	0.00	0.00	0.00	0.00	0.00
120-255	0	0	0.00	0.00	0.00	0.00	0.00
130-255	0	0	0.00	0.00	0.00	0.00	0.00
140-255	0	0	0.00	0.00	0.00	0.00	0.00
150-255	0	0	0.00	0.00	0.00	0.00	0.00
160-255	0	0	0.00	0.00	0.00	0.00	0.00
170-255	0	0	0.00	0.00	0.00	0.00	0.00
180-255	0	0	0.00	0.00	0.00	0.00	0.00
190-255	0	0	0.00	0.00	0.00	0.00	0.00
200-255	0	0	0.00	0.00	0.00	0.00	0.00
210-255	0	0	0.00	0.00	0.00	0.00	0.00
220-255	0	0	0.00	0.00	0.00	0.00	0.00
230-255	0	0	0.00	0.00	0.00	0.00	0.00
240-255	0	0	0.00	0.00	0.00	0.00	0.00
250-255	0	0	0.00	0.00	0.00	0.00	0.00

Table 4b

(continuation, to the right, of Table 4a)

THRESHOLD RANGE	TOTAL PERIMETER POINTS	AVERAGE PERIMETER POINTS PER OBJECT	TOTAL PERIMETER	AVERAGE PERIMETER PER OBJECT	AVERAGE HORIZONTAL FERET DIAMETER	AVERAGE VERTICAL FERET DIAMETER	CUMULATIVE 100	AVERAGE 100 PER OBJECT	AVERAGE CUMULATIVE 100 PER OBJECT
0-255	20	20.00	20.00	20.00	4.00	4.00	700	700.00	30.00
10-255	23	11.50	28.00	14.00	3.30	3.30	700	350.00	40.00
20-255	0	0.00	0.00	0.00	2.50	2.00	420	155.00	37.50
30-255	0	0.00	0.00	0.00	2.23	1.73	480	120.00	40.00
40-255	0	0.00	0.00	0.00	2.00	2.00	300	150.00	30.00
50-255	0	0.00	0.00	0.00	1.00	1.00	100	30.00	30.00
60-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
70-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
80-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
90-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
100-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
110-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
120-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
130-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
140-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
150-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
160-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
170-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
180-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
190-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
200-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
210-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
220-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
230-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
240-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00
250-255	0	0.00	0.00	0.00	0.00	0.00	0	0.00	0.00

C. Literature Review

Emphasis has been placed on generating and updating a cumulative bibliography pertaining to this project (page 27). For literature to be included in the bibliography, it must first satisfy one or more of the following criteria:

1. Alternative and comparable approaches to autoradiography, such as high-speed autoradiography and scintillation autoradiography.
2. Physical modelling, such as the energy transfer between a radionuclide and scintillant, and single photon detection and capture.
3. Specimen preparation techniques, such as animal tissue and glass-slide preparation, availability and assortment of radionuclides and scintillants, and techniques of scintillation counting.
4. Video intensified microscopy, such as photodetectors and their available technology (photomultiplier tubes, charge-coupled devices, microchannel plates), signal processing (noise analysis and filtering), image processing and enhancement, and the analysis, display and communications of data.

IV. DISCUSSION

The Microstructural Analysis Routine (MSAR) was designed and implemented to address a set of experimental needs for this work. Specifications required that program execution time be minimized without additional computing hardware, and that two types of output files be created for each specified threshold range.

Upon program execution, one file contains an array of the uniquely labeled microstructures and, in addition, any information concerning the creation of the file or necessary for image display. Another file contains the calculated statistics for each unique microstructure. A third file is generated which gives the total and average microstructure statistics for the threshold ranges of a particular image. Additional specifications pertaining to the analyses performed by MSAR were also stipulated:


1. Read Image Data into VAX 11/780 Memory
2. Define Region of Interest (ROI)
3. Select Threshold
4. Determine Adjacency Criteria
5. Calculate Microstructure Statistics

These are described in detail in Section III. B.

A test image containing various patterns, created by varying pixel intensity levels within a matrix, was produced to examine the accuracy of the Microstructural Analysis Routine. These patterns simulated microstructures with edge points on the matrix boundary, with diagonal connections, and with holes. Results indicate that microstructures are being detected and analyzed as designed.

Using MSAR, the linearity of the system response as a function of integration time was assessed (Figure 4). In this plot the pixel count above threshold is shown as a function of the integration time. The ratio of photoevent area to full image field is given as a percentage and is printed adjacent to each data point. There is a linearity of the relationship for unsaturated images (i.e., when less than 44.5% of the image field is affected) indicating that, for a fixed source, the integration time is proportional to the cumulative area of the detected photoevents, until saturation effects alter this correlation due to overlapping events. Saturation begins to affect this linear relationship once the photoevent density reaches a concentration greater than 44.5% of the image field (Figure 4). At this value, continuing accumulation results in a progressively higher rate of overlapping points. All points are counted only once, hence leading to an underestimation of photoevents after the saturation point is reached.

FIGURE 4. Plot of Image Area above Threshold versus Signal Integration Time



V. CONCLUSIONS

During this second year, final system evaluations of The Hamamatsu C1966-20 VIM Photonic Microscopy System and the Gould IP8400 Programmable Image Display Unit were carried out. Upon receiving official approval, purchase of the items was ordered.

A Microstructural Analysis Routine (MSAR) was developed to analyze the extremely photon-limited events recorded by the System under development. A program (Particlet) obtained from the United States Army Medical Research Institute of Chemical Defense (USAMRICD) was used as a fundamental concept in successfully accomplishing this task.

Two peer-reviewed scientific manuscripts reporting upon significant aspects of the work to date were prepared and presented, accepted for publication, or published during this reporting period:

"Detection and Localization of Quantum-Limited Events from Radionuclide Labeled Material by Computer Enhanced Video Microscopy." RK Traub, NJ Pressman, JK Frost, PK Gupta, RL Showers, GW Gill, DL Cook and JK Frost, Jr. In the Optical Society of America's Topical Meeting on Quantum Limited Imaging and Image Processing Technical Digest, 86(6):63-66, 1986.

"A Quantitative Method for the Detection and Localization of Quantum-Limited Events From Radionuclides in Cells and Tissue Sections by Computer Enhanced Video Microscopy." NJ Pressman, JK Frost, PK Gupta, RL Showers, GW Gill, DL Cook, JK Frost, Jr and RK Traub. Presented at the International Academy of Cytology's Fifth International Conference on Automation of Diagnostic Cytology and Histology, Brussels, Belgium, May 30-June 1, 1986. *Analyt Quant Cytol Histol*, (accepted for publication).

BIBLIOGRAPHY

Allen RD, Allen NS: Video-enhanced microscopy with a computer frame memory. *J Microsc* 129(Pt 1):3-17, 1983.

Allen RD, Allen NS, Travis JL: Video-enhanced contrast, differential interference contrast (AVEC-DIC) microscopy: a new method capable of analyzing microtubule-related motility in the reticulopodial network of *Allogromia laticollaris*. *Cell Motil* 1:291-302, 1981.

Allen RD, Metuzals J, Tasaki I, Brady ST, Gilbert SP: Fast axonal transport in squid giant axon. *Science* 218:1127-1129, 1982.

Arndt-Jovin DJ, Robert-Nicoud M, Kaufman SJ, Jovin TM: Fluorescence digital imaging microscopy in cell biology. *Science* 230:247-256, 1985.

Arnold BA, Schiebe PO: Noise analysis of a digital radiography system. *Radiology* 142:609-613, 1984.

Axelrod D, Thompson NL, Burghardt TP: Total internal reflection fluorescent microscopy. *J Microsc* 129(Pt 1):19-28, 1983.

Barnard RD: *Thermoelectricity in Metals and Alloys*. John Wiley & Sons, New York, 1972.

Barrows GH, Sisken JE, Allegra JC, Grasch SD: Measurement of fluorescence using digital integration of video images. *J Histochem Cytochem* 32:741-746, 1984.

Brenner M: Computer enhancement of low-light microscopic images. *Am Lab* 15(12):30-35, 1983.

Cowden RR: Fluorescent cytochemistry of living cells. *Am Biotechn Lab* 3(6):8-15, 1985.

Cowen T, Burnstock G: Image analysis of catecholamine fluorescence. *Brain Res Bull* 9:81-86, 1982.

Davidson JB, Case AL: Rapid electronic auto-fluorography of labeled macromolecules on two-dimensional gels. *Science* 215:1398-1400, 1982.

Docchio F, Ramponi R, Sacchi CA, Bottioli G, Freitas I: Time-resolved fluorescence microscopy of hematoporphyrin-derivative in cells. *Lasers Surg Med* 2:21-28, 1982.

Durie BGM: Scintillator distribution in high-speed autoradiography. *Science* 195:208, 1977.

Durie BGM, Salmon SE: High speed scintillation autoradiography. *Science* 190:1093-1095, 1975.

Eastman Kodak Company: Kodak Products for Liquid Scintillation Counting (Dataservice Publication JJ-59), Rochester, 1980.

Feinberg G: Light. *Sci Am* 219(3):50-59, 1968.

Fischer HA: High resolution scintiautoradiography. *J Histochem Cytochem* 27:1527-1528, 1979.

Forman DS, Turriff DE: Video intensification microscopy (VIM) as an aid in routine fluorescence microscopy. *Histochemistry* 71:203-208, 1981.

Fox BW: Techniques of Sample Preparation for Liquid Scintillation Counting — Part I. In *Laboratory Techniques in Biochemistry and Molecular Biology*, Vol. 5, Work TS, Work E (Eds), American Elsevier, New York, 1976, pp 1-333.

Freitas MI, Giordano PA, Bottiroli G: Improvement in microscope photometry by voltage to frequency conversion: analogue measurement and digital processing.

J Microsc 124:211-218, 1979.

Frost JK, Pressman NJ, Showers RL, Gill GW, Albright CD, Frost JK Jr: Characterization of cellular activity by spectral information. *Anal Quant Cytol* 4:146, 1982 (abstract).

Gibson JAB: Liquid Scintillation Counting of Novel Radionuclides. In *Liquid Scintillation Science and Technology*. Noujaim AA, Ediss C, Weibe LI (Eds). Academic Press, New York, 1976, pp 153-166.

Goldgefter L, Toder V: Scintillator distribution in high-speed autoradiography. *Science* 195:208, 1977.

Goldstein DJ: Relation of effective thickness and refractive index to permeability of tissue components in fixed sections. *J R Microsc Soc* 84:43-45, 1965.

Gottlinger C, Meyer KL, Weichel W, Muller W, Raftery B, Radbruch A: Cell-cooling in flow cytometry by Pettier elements. *Cytometry* 7:295-297, 1986.

Grant LD, Stumpf WE: Combined autoradiography and formaldehyde-induced fluorescence methods for localization of radioactively labeled substances in relation to monoamine neurons. *J Histochem Cytochem* 29:175-180, 1981.

Hamamatsu Photonics KK: Characteristics and Applications of Microchannel Plates. Technical Manual RES-0795, Hamamatsu City, Japan, 1982.

Haselmann H: Light microscopy 1981: the state of the art. *J Microsc* 129:1-2, 1983.

Hayakawa T: Ultra-high sensitivity VIM camera system. *Seibutsu Butsuri (Biophysics)* 24:173-178, 1984.

Hayakawa T, Kinoshita K, Miyaki S, Fujiwake H, Ohsuka S: Ultra-low-light level camera for photon counting imaging. *Photochem Photobiol* 41:1-3, 1985.

Hough PVC, McKinney WR, Ledbetter MC, Pollack RE, Moos HW: Identification of biological molecules *in situ* at high resolution via the fluorescence excited by a scanning electron beam. *Proc Natl Acad Sci USA* 73:317-321, 1976.

Howard PL: Basic Liquid Scintillation Counting. American Society of Clinical Pathologists, Chicago, 1976.

Huang, HK, Meyers GL, Martin RK, Albright JP: Digital processing of microradiographic and fluorochromic images from bone biopsy. *Comput Biol Med* 13:27-47, 1983.

Inaba H, Yamagishi A, Takyu C: Development of an ultra-high sensitive photon counting system and its application to biomedical measurements. *Opt Las Eng* 3:125-130, 1982.

Ishida M, Doi K, Loo LN, Metz CE, Lehr JL: Digital image processing: effect on detectability of simulated low-contrast radiographic patterns. *Radiology* 150:569-575, 1984.

King RE: Digital image processing in radioisotope scanning. *IEEE Trans Biomed Eng* 21:414-416, 1974.

Kobayashi Y: LSC Applications Notes #31-50, New England Nuclear, Boston, 1980.

Kobayashi Y, Harris WG: LSC Applications Notes #1-30, New England Nuclear, Boston, 1983.

Kobayashi Y, Maudsley DV: Biological Applications of Liquid Scintillation Counting. Academic Press, New York, 1974.

Kopin IJ, Palkovits M, Kobayashi RM, Jacobowitz DM: Quantitative relationship of catecholamine content and histofluorescence in brain of rats. *Brain Res* 80:229-235, 1974.

Kopriwa BM: Examination of various methods for the intensification of radioautographs by scintillators. *J Histochem Cytochem* 27:1524-1526, 1979.

Kopriwa BM: Quantitative investigation of scintillator intensification for light and electron microscope radioautography. *Histochemistry* 68:265-279, 1980.

Kristinn J, Blouke M: Charge-coupled devices in astronomy. *Sci Am* 247(4):67-74, 1982.

Langager JM, Howard GA, Baylink DJ: An improved technique for rapid autoradiography of cells and tissue sections. *Histochemistry* 75:523-531, 1982.

Mangravite JA, Gallis D, Foery R: The recovery of organic solvents from liquid scintillation waste. *Am Lab* 15(7):24-34, 1983.

Mankovich NJ, Trambert MA, Huang HK: The architecture of a prototype multimodality image viewing station. In *Proceedings of the Conference on Medical Image Production, Processing, Display, and Archiving*. San Diego, SPIE, 1984.

Matsuura S, Umebayashi S, Okuyama C, Oba K: Characteristics of the newly developed MCP and its assembly. *IEEE Trans Nucl Sci* NS32:350-354, 1985.

Miller M: Specialized video operations. *Video Systems* 10(Sep):32-47, 1984.

Miller TR, Sampathkumaran KS, King MA: Rapid digital filtering. *J Nucl Med* 24:625-628, 1983.

McGregor GN, Kapitza HG, Jacobson KA: Laser-based fluorescence microscopy of living cells. *Laser Focus* 20(Nov):85-93, 1984.

New England Nuclear: *Liquid Scintillation Counting — Products and Procedures*, Boston, 1982.

O'Kane DJ, Palevitz BA: Video microscopy for low-light objects. *J Cell Biol* 83:304A, 1979.

Panayi GS, Neill WA: Scintillation autoradiography—a rapid technique. *J Immunol Methods* 2:115-117, 1972.

Peng CT, Horrocks DL, Alpen EL (Eds): *Liquid Scintillation Counting — Recent Applications and Development*, Vols. I and II. Academic Press, New York, 1980.

Pressman NJ: Implementation of an automated high-resolution audio-video teleconferencing system for biomedical image communications, processing and analysis. *Comput Biol Med* (in preparation).

Pressman NJ, Frost JK: BIDAS: A microscope imaging, communications, processing and analysis system for cytopathology. In *Proceedings of the IEEE 1984 International Symposium on Medical Images and Icons*. Silver Spring, Maryland, IEEE Computer Society Press, 1984, pp 372-379.

Pressman NJ, Frost JK: Design of a high-resolution microcomputer-controlled audio-video teleconferencing system for medical education, research and patient care. *Diagn Cytopathol* (in preparation).

Pressman NJ, Frost JK, Gupta PK, Showers RL, Gill GW, Cook DL, Frost, JK Jr, Traub RK: The Development of a System for Enhanced Radionuclide Detection, Annual Report, August 1, 1984-July 31, 1985, Quantitative

Cytopathology Laboratories, Cytopathology Division, Department of Pathology, The Johns Hopkins University School of Medicine, Baltimore, Maryland. Supported by U. S. Army Medical Research and Development Command Contract No. DAMD17-84-C-4222.

Pressman NJ, Frost JK, Gupta PK, Showers RL, Gill GW, Cook DL, Frost JK Jr, Traub RK: A quantitative method for the detection and localization of quantum limited events from radionuclides in cells and tissue sections by computer enhanced video microscopy. Accepted for presentation at the Fifth International Conference on Automation of Diagnostic Cytology and Histology, Hilton International Hotel, Brussels, Belgium, May 30-June 1, 1986. *Anal Quant Cytol Histol* (accepted for publication in 1987).

Pressman NJ, Frost JK, Showers RL, Frost JK Jr, Gill GW, Albright CD: Measurement of spectral energy absorption in cytologic preparations using scanning transmission optical microscopy. *Anal Quant Cytol* 4:155, 1982 (abstract).

Pressman NJ, Frost JK, Showers RL, Frost JK Jr, Gill GW, Cook DL, Tyrer HW: The Johns Hopkins scanning transmission optical microscope (STOM): I. overview of design, implementation and performance for automated cytopathology. *Rev Sci Instrum* (in preparation).

Pressman NJ, Frost JK Jr, Showers RL, Frost JK: Design of an interactive color-image cytometry system for data acquisition, communications, processing and analysis. *Pathol Res Pract* (in preparation).

Pressman NJ, Frost JK Jr, Showers RL and Frost JK: Development and implementation of an interactive color-image cytometry system. *Anal Quant Cytol Histol* (in preparation).

Pressman NJ, Showers RL, Gill GW, Frost JK Jr, Cook DL, Gupta PK, Frost JK, Traub RK: Computer enhanced video microscopy for quantitative analyses of low concentrations of radionuclide labeled compounds in cells and tissues. *J Opt Soc Am* (in preparation).

Pressman NJ, Tyrer HW, Frost JK: An analytical cytology video-audio network for teleconferencing and computing. *Cytometry* (in preparation).

Preston K Jr: High-resolution image analysis. *J Histochem Cytochem* 34:67-74, 1986.

Prudhomme D: A cooling device for flow Cytometric Systems. *Stain Technol* 53:300-301, 1978.

Przybylski RJ: Scintillation radioautography: a new technique designed to augment silver grain number in radioautographs. *J Cell Biol* 43:108a, 1969.

Reynolds GT: Image intensification applied to microscope systems. In *Advances in Optical and Electron Microscopy*, Vol. 2. Academic Press, New York, 1968, pp 1-40.

Reynolds, GT: Applications of image intensification to low level fluorescence studies of living cells. *Microsc Acta* 83:55-62, 1980.

Reynolds GT, Taylor DL: Image intensification applied to light microscopy. *BioScience* 30:586-592, 1980.

Rich ES Jr, Wampler JE: A flexible, computer-controlled video microscope capable of quantitative, spatial, temporal, and spectral measurements. *Clin Chem* 27:1558-1568, 1981.

Riva RD: Exfoliative cytology—a review of literature and presentation of a fluorescent technique. *J Oral Med* 32(4):98-103, 1977.

Robinson AL: Demonstrating single photon interference. *Science* 231:671-672, 1986.

Rosenquist TH: On the use of quantitation in histochemistry: a review of the literature. *Histochem J* 8:205-208, 1976.

Roth LJ, Diab IM, Watanabe M, Dinerstein RJ: A correlative radioautographic, fluorescent, and histochemical technique for cytopharmacology. *Mol Pharmacol* 10:986-998, 1974.

Sedlacek W, Gundlach H, Ax W: The use of television cameras equipped with an image intensifier in immunofluorescence microscopy. *Behring Inst Mitt* 59:64-70, 1976.

Siegel JI, Shumsky E: Image analysis and quantitative microscopy. *Photomethods* 29(2):34-36, 1986.

Silage DA, Gil J: Digital reversal and enhancement of negatives in video-based interactive morphometry. *Lab Invest* 51(1):112-116, 1984.

Sisken JE, Barrows GH, Grasch SD: The study of fluorescent probes by quantitative video intensification microscopy (QVIM). *J Histochem Cytochem* 34:61-66, 1986.

Smith DA, Shelley SA, Balis JU: Reduced cell proliferation in fetal lung after maternal administration of pilocarpine: a scintillation autoradiographic study. *Am J Anat* 155:131-137, 1979.

Stanulis BM, Sheldon S, Grove GL, Cristofalo VJ: Scintillation fluid shortens exposure times in autoradiography. *J Histochem Cytochem* 27:1303-1307, 1979.

Sullivan JA: Computer-enhanced video microscopy. *Photomethods* 29(2):31-33, 48, 1986.

Tait WH: Unsharpness and contrast in digitized images. *Phys Med Biol* 29:1177-1186, 1982.

Tanasugarn L, McNeil P, Reynolds GT, Taylor DL: Microspectrofluorometry by digital image processing: measurement of cytoplasmic pH. *J Cell Biol* 98:717-724, 1984.

Taylor, DL: Fluorescence Microscopy. Course notes from Analytical and Quantitative Light Microscopy in Biology, Medicine, and Material Sciences, private communication. Seminar directed by Shinya Inoue, Woods Hole, Massachusetts, March, 1984 (unpublished).

Taylor DL, Wang YL: Fluorescently labeled molecules as probes of the structure and function of living cells. *Nature* 284:405-410, 1980.

Trager W, Stanley HS, Allen RD, Allen NS: Knobs on the surface of erythrocytes infected with plasmodium falciparum: visualization by video-enhanced, differential interference contrast microscopy. *J Parasitol* 68:332-333, 1982.

Traub K: In vivo distribution of ^{14}C radiolabeled Soman ((3,3-dimethyl-2-butoxy)-methylphosphorylfluoride) in the central nervous system of the rat. *Neurosci Lett*, NSL-03542, 60:219-225, 1985.

Traub RK, Pressman NJ, Frost JK, Gupta PK, Showers RL, Gill GW, Cook DL, Frost JK Jr: Detection and localization of quantum limited events from radionuclide labeled materials by computer enhanced video microscopy. In the Optical Society of America's Topical Meeting on Quantum Limited Imaging and Image Processing Technical Digest, 86(6):63-66, 1986.

Travis JL, Allen RD: Studies on the motility of the foraminifera I. ultrastructure of the reticulopodial network of *Allogromia laticollaris* (Arnold). *J Cell Biol* 90:211-221, 1981.

Van der Ploeg M, Bosman FT, Van Duijn P: Pre-exposure of films in fluorescence photography. *Histochem J* 8:201-204, 1976.

Van der Ploeg M, Van Duijn P: Reflection versus fluorescence. *Histochemistry* 62:227-232, 1979.

Walter RJ, Berns MW: Computer-enhanced video microscopy: digitally processed microscope images can be produced in real time. *Proc Natl Acad Sci USA* 78:6927-6931, 1981.

Walton AJ, Templer RH, Reynolds GT: Depth of field measurements relevant to single photon detection using image-intensified microscopy. *Optica Acta* 32:1349-1360, 1986.

Wick RA: Grabbing images at very low light levels. *Photonics Spectra*, 19(5): 133-136, 1985.

Wied GL, Bibbo M, Bartels PH: Computer analysis of microscopic images: application in cytopathology. *Pathol Annu* 16(Pt 1):367-409, 1981.

Williams MA: Autoradiography: its methodology at the present time. *J Microsc* 128(Pt 1):79-94, 1982.

Woodcock CLF, D'Amico-Martel A, McInnis CJ, Annunziato AT: How effective is 'high-speed' autoradiography? *J Microsc* 117(Pt 3):417-423, 1979.

Young IT, Balasubramanian, Dunbar DL, Peverini RL, Bishop RP. SSAM: Solid-state automated microscope. *IEEE Trans Biomed Eng* 29:70-81, 1982.

Young IT, Paskowitz L: Localization of cellular structures. *IEEE Trans Biomed Eng* 22:35-40, 1975.

GLOSSARY


OF SELECTED KEY TERMS, UNITS, AND ABBREVIATIONS

AChE	-Acetylcholinesterase
attogram	- 10^{-18} gram
ASA	- American Standards Association, speed of photographic film
BIDAS	- The Johns Hopkins Biomedical Image Data Analysis System
Ci	- Curie = 3.7×10^{10} disintegrations per second
^{14}C	- Radioactive Carbon isotope with atomic mass of 14
C1966-20	- Hamamatsu VIM image processing system
CFIAC	- Cytotechnologist Fellow, International Academy of Cytology
cpm	- Counts per minute
cps	- Counts per second
DAGE-MTI	- Dage-MTI, Inc., Michigan City, IN 46360
DAMD17-84-C-4222	- USAMRICD Contract for "The Development of a System for Enhanced Radionuclide Detectors" with JHU
DEC	- Digital Equipment Corporation, Northboro, MA 01532
dps	- Disintegrations per second
EM	- Electron Microscopy
femtogram	- 10^{-15} gram
feret diameter	- The maximum diameter of a structure in a given direction
FIAC	- Fellow, International Academy of Cytology
FITC	- Fluorescein Isothiocyanate
flux	- Rate of transfer of particles or energy per unit area
GOULD	- Gould Inc., DeAnza Imaging & Graphics Division, San Jose, CA 95131
h	- Planck's constant = 6.62363×10^{-34} joule-second
^3H	- Isotope of Hydrogen with atomic mass of 3
IAC	- International Academy of Cytology
IOD	- Integrated Optical Density
JHMI	- The Johns Hopkins Medical Institutions
JHU	- The Johns Hopkins University
mCi	- Millicurie = 10^{-3} curie
MCP	- Microchannel Plate

MIAC	- Member, International Academy of Cytology
MSAR	- The Microstructural Analysis Routine
OD _{max}	- Maximum Optical Density
OSA	- Optical Society of America
Particlet	- Particle counting and sizing program developed by USAMRICD
P ² /A	- Perimeter squared divided by Area
PI	- Principal Investigator
pixel	- Picture element
PMI	- Photonic Microscopy, Inc., Oak Brook, IL 60521
PMIAC	- Professional Member, International Academy of Cytology
RFP	- Request for Proposals
RGB	- Red/Green/Blue
ROI	- Region of Interest
slice-value	- Boundary value for a threshold range
specific activity	- Activity (i.e., disintegrations per second) in (micro-, milli-) curies of radioactive species per unit mass or concentration of non-radioactive material (e.g., 56 mc C ¹⁴ /mM experimental material)
S/N	- Signal to Noise ratio
³⁵ S	- Isotope of Sulfur with atomic mass of 35
USAMRICD	- U. S. Army Medical Research Institute for Chemical Defense
USAMRDC	- U. S. Army Medical Research and Development Command
VAX-11/780	- DEC 32-bit computer
VENUS	- Zeiss low-light-level video camera
VIM	- Video Intensified Microscopy
ZEISS	- Carl Zeiss, Inc., Thornwood, New York 10594


APPENDICES

Appendix A. MSAR Flowcharts




MAIN PROGRAM: EDC_MSAR

Page 34




SUBROUTINE: NFNAME

-SUPPLY IMAGE FILENAME
AND READ IMAGE INTO
VAX MEMORY.




SUBROUTINE: DEFROI

-ALTER REGION OF INTEREST
(ROI) OF IMAGE THAT THE
CALCULATIONS ARE
PERFORMED WITHIN.




SUBROUTINE: SETTHR

-ALTER THRESHOLD VALUES.
THESE DETERMINE RANGE
OF PIXEL VALUES
CONSIDERED IN ANALYSIS.




SUBROUTINE: ADJCRI

-ALTER RULE GOVERNING
CRITERIA FOR DETERMINING
ADJACENCY OF
NEIGHBORING PIXELS.




SUBROUTINE: EXECUTE

-PERFORM MICROSTRUCTURAL
DETECTION AND ANALYSIS.



SUBROUTINE: CONTIG


-DETERMINE ADJACENCY
OF PIXELS AND
ASSIGN PIXELS TO
DISTINCT
MICROSTRUCTURES.



SUBROUTINE: STATIS


-CALCULATE STATISTICS
ON EACH
MICROSTRUCTURE.

Page 41



SUBROUTINE: ANLZST

-MICROSTRUCTURAL
DATA ANALYSIS.



Appendix B. MSAR Program Code

```

PROGRAM MSAR
C
C... Read in a REAL *4 image from VAX, measure microstructure
C... information (e.g. number, geometric X&Y center, weighted X&Y center,
C... area, perimeter length, perimeter count, maximum horizontal and
C... vertical feret diameters, Integrated Optical Density and average
C... Optical Density).
C
      INCLUDE 'EDC_MSARCOM.FOR'
      BYTE INFTXT(80)
      DATA IOPT/1/
C
C... PRINT RUNNING INFO TO USER
C
      OPEN(UNIT=1, NAME='[EDC.TASKS]EDC_MSAR_TEXT.FOR', TYPE='OLD',
C READONLY)
5      READ(1,10,END=15), NCHAR, (INFTXT(J), J=1, NCHAR)
10     FORMAT(0, <NCHAR>A1)
      TYPE 12, (INFTXT(J), J=1, NCHAR)
12     FORMAT(' ' <NCHAR>A1)
      GOTO 5
15     CLOSE(UNIT=1, DISPOSE='SAVE')
      TYPE 16
16     FORMAT('$TYPE CARRIAGE RETURN TO CONTINUE : ')
      ACCEPT 17, INFTXT(1)
17     FORMAT(A1)
C
C
C... DISPLAY INFO TO USER IN MENU SCREEN.
C
      CALL WRTTMP
      CALL WRTINF
C
C... DISPLAY MENU TO USER AT BOTTOM OF SCREEN
C
20     CALL WRTMEN
C
C... DISPLAY MENU SELECTION CURSOR. WAIT FOR USER TO SELECT OPTION
C
      CALL SELOPT(IOPT)
      GOTO (100, 200, 300, 400, 500, 600), IOPT
      GOTO 20
C
100    CALL NFNAME                                !SET NEW FILENAME
      GOTO 20
200    IF(FNMSET .EQ. 0) GOTO 900                  !CHECK IF FILENAME NOT SPECIFIED
      CALL DEFROI                                  !DEFINE ROI
      GOTO 20
300    CALL ADJCRI                                  !ADJACENCY CRITERIA
      GOTO 20
400    CALL SETTHR                                  !SET THRESHOLDS
      GOTO 20
500    IF(FNMSET .EQ. 0) GOTO 900                  !CHECK IF FILENAME NOT SPECIFIED
      CALL EXECUTE                                !EXECUTE
      GOTO 20
600    CALL CLRSCR                                  !EXIT FROM ENTIRE PROGRAM
      CALL CLRTXT
      STOP
C
C... OPERATION CANNOT BE PERFORMED UNLESS FILE HAS BEEN READ INTO MEMORY

```

```

C
900  CALL CLRLIN(19,23)          !CLEAR LINES 19-23 (MENU LOCATION)
      CALL SCROLL(19,23)        !SET SCROLL REGION FOR 19-23
      TYPE 910
910  FORMAT('/' OPERATION CANNOT BE PERFORMED UNTIL FILE HAS BEEN
C    READ INTO MEMORY'/)
      T=SECNDS(0.)
920  CONTINUE
      IF(SECNDS(0.)-T .LT. 2.0) GOTO 920
      CALL SCROLL(1,24)
      CALL CLRLIN(19,23)
      GOTO 20
C
      END
C
C
C
C

```

```

SUBROUTINE WRTTMP
C
C... WRITE TEMPLATE TO TERMINAL
C
CALL CLRTXT
CALL VTCURP(4,1)
TYPE 100
TYPE 120
TYPE 130
TYPE 140
TYPE 150
TYPE 160
TYPE 170
TYPE 180
TYPE 185
TYPE 190
TYPE 200
100 FORMAT('+FILENAME: '//)
120 FORMAT(' ')
130 FORMAT(' THRESHOLD RANGE
140 FORMAT(' LOWER: UPPER:
C          COLUMNS: ')
150 FORMAT(' THRESHOLD INCREMENT:
C          Y-TOP: ')
160 FORMAT(' LOWER_INCR: UPPER_INCR:
C          Y-BOTTOM: ')
170 FORMAT(' ')
180 FORMAT(' ADJACENCY CRITERIA: ')
185 FORMAT(' PIXEL DIMENSION: MICRONS')
190 FORMAT(' ')
200 FORMAT(' PROGRAM STATUS: ')
C
RETURN
END
C
C
C
C
C

```

```

SUBROUTINE WRITMEN
C
C... WRITE MENU TO TERMINAL
C
CHARACTER*10 OPTION(6)
INTEGER *2 OPTLEN(6),OPTLIN(6),OPTCOL(6)
BYTE BOLD(4),ATTOFF(4),REVERS(4)
DATA OPTION(1) //'FILENAME '// !MENU OPTIONS
DATA OPTION(2) //'DEFINE ROI'//
DATA OPTION(3) //'ADJACENCY '//
DATA OPTION(4) //'THRESHOLD '//
DATA OPTION(5) //'EXECUTE '//
DATA OPTION(6) //'EXIT '//
DATA OPTLEN/8,10,9,9,7,4/ !LENGTH OF MENU OPTIONS
DATA OPTLIN/20,20,20,22,22,22/
DATA OPTCOL/1,20,39,1,20,39/
DATA BOLD/27,91,49,109/ !ESC [ 1 m
DATA ATTOFF/27,91,48,109/ !ESC [ 0 m
DATA REVERS/27,91,55,109/ !ESC [ 7 m
C
TYPE 90, (BOLD(J), J=1, 4)
90 FORMAT('+'A1)
DO 100 IOPT=1,6
CALL VTCURP(OPTLIN(IOPT),OPTCOL(IOPT))
TYPE 95,OPTION(IOPT)(1:10)
95 FORMAT('+'A10)
100 CONTINUE
RETURN
END
C
C
C
C
C

```

```

SUBROUTINE SELOPT(IOPT)
CHARACTER*10 OPTION(6)
INTEGER *2 OPTLEN(6), OPTLIN(6), OPTCOL(6)
BYTE BOLD(4), ATTOFF(4), REVERS(4), BLINK(4)
DATA OPTION(1) /'FILENAME' //          !MENU OPTIONS
DATA OPTION(2) /'DEFINE ROI'//
DATA OPTION(3) /'ADJACENCY'//
DATA OPTION(4) /'THRESHOLD'//
DATA OPTION(5) /'EXECUTE'//
DATA OPTION(6) /'EXIT'//
DATA OPTLEN/8, 10, 9, 9, 7, 4/          !LENGTH OF MENU OPTIONS
DATA OPTLIN/20, 20, 20, 22, 22, 22/
DATA OPTCOL/1, 20, 39, 1, 20, 39/
DATA BOLD/27, 91, 49, 109/             !ESC [ 1 m
DATA ATTOFF/27, 91, 48, 109/           !ESC [ 0 m
DATA REVERS/27, 91, 55, 109/           !ESC [ 7 m
DATA BLINK/27, 91, 53, 109/            !ESC [ 5 m

C
JOPT=IOPT

C
60 TYPE 90, (ATTOFF(J), J=1, 4)          !RETURN 'OLD' OPTION TO NORMAL
TYPE 90, (BOLD(J), J=1, 4)              !HOWEVER, BOLD THE OPTION
CALL VTCURP(OPTLIN(JOPT), OPTCOL(JOPT))
TYPE 95, OPTION(JOPT)(1:10)

C
TYPE 90, (REVERS(J), J=1, 4)            !'REVERSE' 'NEW' OPTION
C
TYPE 90, (BLINK(J), J=1, 4)             !'BLINK' 'NEW' OPTION
90 FORMAT('+ 'A1)
CALL VTCURP(OPTLIN(IOPT), OPTCOL(IOPT))
TYPE 95, OPTION(IOPT)(1:10)
95 FORMAT('+ 'A10)
CALL VTCURP(23, 1)

C
C... WAIT FOR USER TO MOVE MENU CURSOR
C

IKEY=0
TLAST=SECNDS(0.)                        !SECONDS SINCE LAST INPUT
100 I=IESC(IKEY)
IF(IKEY .EQ. 0) THEN
    IF(SECNDS(0.)-TLAST .GT. 3600.) THEN
        IOPT=6
        GOTO 999
    END IF
    GOTO 100
END IF
IF(IKEY .EQ. 13) GOTO 999                !<CR> USER ACCEPTS IOPT
IF(IKEY .LT. 65 .OR. IKEY .GT. 68) GOTO 100
JOPT=IOPT
IF(IKEY .EQ. 65) THEN                    !UP ARROW
    IOPT=IOPT-3
    IF(IOPT .LT. 1) IOPT=IOPT+6
END IF
IF(IKEY .EQ. 66) THEN                    !DOWN ARROW
    IOPT=IOPT+3
    IF(IOPT .GT. 6) IOPT=IOPT-6
END IF
IF(IKEY .EQ. 67) THEN                    !RIGHT ARROW
    IOPT=IOPT+1
    IF(IOPT .GT. 6) IOPT=1
END IF

```

```
IF(IKEY .EQ. 68) THEN          !LEFT ARROW
  IOPT=IOPT-1
  IF(IOPT .LT. 1) IOPT=6
END IF
GOTO 60
```

```
C
999 TYPE 90, (ATTOFF(J), J=1,4)  !RETURN OPTION TO NORMAL
RETURN
END
```

```
C
C
C
C
```

```

SUBROUTINE NFNAME
INCLUDE 'EDC_MSARCOM.FOR'

C
C...  SETUP SOME PARAMETERS
C
      IRROW=512
      IRCOL=512
      IRXTOP=1
      IRYTOP=1
      IRXBOT=512
      IRYBOT=512

C
C...  CLEAR MENU SECTION OF VT SCREEN
C
      CALL CLRLIN(19,23)          !CLEAR LINES 19-23 (MENU LOCATION)
      CALL SCROLL(19,23)         !SET SCROLL REGION FOR 19-23

C
C...  PROMPT USER FOR FILENAME
C
      DO 90 J=1,14
      FNAME(J)=0
      TYPE 100
100    FORMAT('ENTER FILENAME AND EXTENSION : ')
      ACCEPT 120,NCHAR,(FNAME(J),J=1,NCHAR)
120    FORMAT(0,13A1)
      IF(NCHAR.EQ. 0) THEN
          CALL SCROLL(1,24)
          CALL CLRLIN(19,23)
          GOTO 999
      END IF

C
C...  GET PIXEL DIMENSIONS FROM USER
C
      TYPE 320
320    FORMAT(' ENTER PHYSICAL DIMENSION OF PIXEL IN MICRONS
C (ie 0.25)')
      TYPE 330
330    FORMAT('TYPE <CR> IF UNKNOWN : ')
      ACCEPT 340,NCHAR,SIZPIX
340    FORMAT(0,F15.5)
      IPIXFG=0          !SET PIXEL FLAG=0
      IF(NCHAR.EQ. 0) THEN
          IPIXFG=1      !SET PIXEL FLAG SINCE USER DIDN'T KNOW SIZE
          SIZPIX=1.0
      END IF

C
C...  OPEN INPUT FILE
C
      OPEN(UNIT=3,NAME=FNAME,TYPE='OLD',READONLY,FORM='UNFORMATTED',
CERR=900)

C
C...  READ IMAGE INTO MEMORY
C
      TYPE 350
350    FORMAT('///// READING IMAGE FILE INTO MEMORY')
      READ (3) TEMP
      CLOSE(UNIT=3,DISPOSE='SAVE')
      NCOLS=512
      NROWS=512
      DO 400 J=1,512          !TRANSFER TO I*2 ARRAY

```

```

DO 400 K=1,512
IMAGE(K,J)=IINT(TEMP(K,J))
C
C... WRITE NEW INFORMATION TO TERMINAL
C
CALL SCROLL(1,24)
CALL WRTINF          !WRITE INFO TO SCREEN
CALL CLRLIN(19,23)
FNMSET=1             !SET FLAG TO INDICATE FILE READ SUCCESSFULLY
C
999 RETURN
C
C
900 TYPE 910
910 FORMAT(' ERROR IN OPENING FILE. PLEASE REENTER')
GO TO 80
END
C
C
C
C
C

```

```

SUBROUTINE DEFROI
INCLUDE 'EDC_MSARCOM.FOR'

C
C... CLEAR MENU SECTION OF VT SCREEN
C
CALL CLRRLIN(18,23)          !CLEAR LINES 18-23 (MENU LOCATION)
CALL SCROLL(18,23)          !SET SCROLL REGION FOR 18-23

C
C... KEYPAD ENTRY
C
200 TYPE 220
220 FORMAT('ENTER X,Y COORDINATES OF TOP LEFT CORNER OF ROI : ')
ACCEPT 240,NCHAR,11,12
240 FORMAT(G,214)
IF(NCHAR.EQ. 0) GOTO 999
IF(11.LT. 0 .OR. 12.LT. 0) GOTO 300
IF(11.GT. NCOLS .OR. 12.GT. NROWS) GOTO 300

C
250 TYPE 260
260 FORMAT('ENTER X,Y COORDINATES OF BOTTOM RIGHT CORNER OF ROI : ')
ACCEPT 240,NCHAR,13,14
IF(NCHAR.EQ. 0) GOTO 999
IF(13.LT. 11 .OR. 14.LT. 12) GOTO 320
IF(13.GT. NCOLS .OR. 14.GT. NROWS) GOTO 300

C
IRXTOP=11
IRYTOP=12
IRXBOT=13
IRYBOT=14
IRROW=IRYBOT-IRYTOP+1
IRCOL=IRXBOT-IRXTOP+1
GOTO 999

C
C... ERRORS IN ENTERING ROI INFORMATION
C
300 TYPE 310
310 FORMAT(' BOUNDARY DEFINITION ARE OUTSIDE OF LEGAL RANGE')
GOTO 200
320 TYPE 330
330 FORMAT(' BOTTOM RIGHT MUST BE BELOW AND TO THE RIGHT OF TOP LEFT')
GOTO 250

C
999 CALL SCROLL(1,24)
CALL WRITNF          !WRITE INFO TO SCREEN
CALL CLRRLIN(18,23)

C
RETURN
END

C
C
C
C
C

```

```

SUBROUTINE SETTHR
INCLUDE 'EDC_MSARCOM.FOR'
BYTE ACCEPT

C
C... CLEAR MENU SECTION OF VT SCREEN
C
CALL CLRLIN(19,23)          !CLEAR LINES 19-23 (MENU LOCATION)
CALL SCROLL(19,23)         !SET SCROLL REGION FOR 19-23

C
C... GET THRESHOLD VALUES FROM USER
C
90  TYPE 100
100 FORMAT('ENTER INITIAL LOWER,UPPER THRESHOLD VALUES : ')
    ACCEPT 120,NCHAR,11,12
120 FORMAT(0,214)
    IF(NCHAR.EQ. 0) GOTO 999
    IF(11.LT. 0.OR. 12.LT. 11) GOTO 800
    IF(12.GT. 255.OR. 11.GT. 12) GOTO 800

C
    TYPE 160
160 FORMAT('ENABLE THRESHOLD INCREMENTATION ? <N/Y> : ')
    ACCEPT 170,NCHAR,ACCEPT
170 FORMAT(0,A1)
    IF(ACCEPT.EQ. 'Y') GOTO 190

C
C... USER DID NOT WANT THRESHOLD INCREMENTATION. SET VALUES AND RETURN
C
    ITLOW=11
    ITUPR=12
    ITINCR=0
    ITLINC=0
    ITUINC=0
    GOTO 999

C
C... GET INCREMENTATION VALUES FROM USER
C
190  TYPE 200
200  FORMAT('ENTER LOWER,UPPER THRESHOLD INCREMENT VALUES : ')
    ACCEPT 120,NCHAR,13,14
    IF(NCHAR.EQ. 0) GOTO 999
    IF(13.GT. 255.OR. 14.GT. 255) GOTO 820
    IF(13.EQ. 0.AND. 14.EQ. 0) GOTO 840

C
C... SET VALUES AND RETURN
C
    ITLOW=11
    ITUPR=12
    ITINCR=1
    ITLINC=13
    ITUINC=14
    GOTO 999

C
800  TYPE 810
810  FORMAT('// INVALID THRESHOLD VALUES. PLEASE RE-ENTER//')
    GOTO 90

C
820  TYPE 830
830  FORMAT('// THRESHOLD INCREMENT VALUE CANNOT EXCEED 255. PLEASE
C RE-ENTER//')
    GOTO 190

```

```
C
840 TYPE B50
850 FORMAT('// INVALID INCREMENT VALUES. PLEASE RE-ENTER')
      GOTO 190
C
999  CALL SCROLL(1,24)
      CALL WRTINF          !WRITE INFO TO SCREEN
      CALL CLRLIN(18,23)
C
      RETURN
      END
C
C
C
C
C
```

```

SUBROUTINE ADJCRI
INCLUDE 'EDC_MSARCOM.FOR'
C
C... CLEAR MENU SECTION OF VT SCREEN
C
      CALL CLRLIN(19,23)          !CLEAR LINES 19-23 (MENU LOCATION)
      CALL SCROLL(19,23)         !SET SCROLL REGION FOR 19-23
C
C... GET USER INPUT AND CHANGE ADJACENCY CRITERIA
C
      TYPE 100
100   FORMAT('/' SELECT CRITERIA BY WHICH POINTS ARE CONSIDERED
C CONTIGUOUS/' ' 1-ORTHOGONAL POINTS ONLY/' ' 2-ORTHOGONAL AND
C DIAGONAL POINTS')
      TYPE 120
120   FORMAT('ENTER SELECTION : ')
      ACCEPT 140,NCHAR,ISEL
140   FORMAT(0,12)
      IF(NCHAR.EQ. 0) GOTO 999
      IF(ISEL.LT. 1.OR. ISEL.GT. 2) GOTO 90
      IF(ISEL.EQ. 1) IADJCT=0
      IF(ISEL.EQ. 2) IADJCT=1
C
999   CALL SCROLL(1,24)
      CALL WRTINF                !WRITE INFO TO SCREEN
      CALL CLRLIN(19,23)
C
      RETURN
      END
C
C
C
C
C

```

```

SUBROUTINE EXECUTE
INCLUDE 'EDC_MSARCOM.FOR'

C
C... SAVE INITIAL THRESHOLD VALUES
C
JTLOW=ITLOW
JTUPR=ITUPR
NUMAVG=0

C
C... CHANGE PROGRAM STATUS TO 'EXECUTING'
C
STATUS=2
CALL WRTINF

C
C... OPEN OUTPUT FILE TO CONTAIN 'AVERAGE' INFORMATION OF MICROSTRUCTURES
C
OPEN(UNIT=2, NAME='MICROS.DAT', TYPE='UNKNOWN', ACCESS='APPEND')
DO 10 J=1, 14
IF(FNAME(J).EQ. '. ') GOTO 15
CONTINUE
NCHAR=14
GOTO 20
15
NCHAR=J*3
20
WRITE(2, 25), (FNAME(J), J=1, NCHAR), IRXTOP, IRYTOP, IRXBOT, IRYBOT
25
FORMAT(// 'FILENAME = ' <NCHAR>A1// ' REGION OF INTEREST -
CTOP LEFT (X='13', Y='13'), BOTTOM RIGHT (X='13', Y='13')//)
WRITE(2, 30) !PUT OUT HEADER OF FIRST SECTION
30
FORMAT('
C
C'
C
C' AVERAGE AVERAGE AVERAGE AVERAGE//
C' THRESHOLD TOTAL WEIGHTED WEIGHTED//
C' X- Y- X- Y-//
C' RANGE AREA (MICROSTRUCTURES) AREA
C' CENTER CENTER CENTER CENTER//')

C
C
C... CLEAR MENU SECTION OF VT SCREEN
C
CALL CLRRLIN(19, 23) !CLEAR LINES 19-23 (MENU LOCATION)
CALL SCROLL(19, 23) !SET SCROLL REGION FOR 19-23

C
C... GO THROUGH THE INPUT IMAGE FILE IN ARRAY 'IMAGE', THRESHOLD USING
C... THE USER-SUPPLIED VALUES, AND PLACE THE BINARY RESULT IN 'THRESH'.
C... ONLY BE CONCERNED WITH THE REGION OF INTEREST.
C
50
DO 100 I=IRYTOP, IRYBOT
TYPE 99, 1
99
FORMAT(' THRESHING LOC '13)
DO 100 J=IRXTOP, IRXBOT
IF(IMAGE(J, I) .GE. ITLOW .AND. IMAGE(J, I) .LE. ITUPR) THEN
THRESH(J, I)=1
C... *** NOTE - NEXT LINE OF CODE IS FOR VERTICAL FILL OF IMAGE.
IF(IMAGE(J, I-1) .LT. ITLOW .OR. IMAGE(J, I-1)
C.GT. ITUPR) THRESH(J, I-1)=-1
ELSE
THRESH(J, I)=0
END IF
100
CONTINUE
C

```

```

C...   GO BACK INTO IMAGE AND SEARCH FOR -1 FLAG.  IF FOUND, PUT THE VALUE
C...   OF THE PIXEL DIRECTLY VERTICAL TO THE PRESENT PIXEL INTO PRESENT
C...   LOCATION.
C
      DO 2100 I=IRYTOP,IRYBOT
      TYPE 99,I
      DO 2100 J=IRXTOP,IRXBOT
      IF(THRESH(J,I) .EQ. -1) THEN
        THRESH(J,I)=1
        IMAGE(J,I)=IMAGE(J,I+1)
      END IF
2100   CONTINUE
C
C...   CALL SUBROUTINE 'CONTIG' TO CALCULATE THE CONTIGUOUS STRUCTURES
C...   IN THE USER-DEFINED ROI AND USING THE USER-SPECIFIED ADJACENCY
C...   CRITERIA.  RESULT OF THIS SUBROUTINE WILL BE UNIQUELY LABEL
C...   MICROSTRUCTURES IN THE ARRAY 'THRESH'.
C
      CALL CONTIG
      TYPE 101,UNIQUE
101     FORMAT(' UNIQUE NUMBER = 'I7)
      T=SECNDS(0.)
102     CONTINUE
      IF(SECNDS(0.)-T .LT. 2.0) GOTO 102
C
C...   CALL SUBROUTINE 'STATIS' TO CALCULATE THE STATISTICS ON EACH
C...   OF THE UNIQUE MICROSTRUCTURES IN THE ARRAY 'THRESH'.  WRITE
C...   OUTPUT TO OUTPUT FILE
C
      CALL STATIS(NUMAVG)
C
C...   WAS THRESHOLDING INCREMENTATION ENABLED?
C
      IF(ITINCR .EQ. 0) GOTO 900
C
C...   INCREMENT THRESHOLDS
C
      ITLOW=ITLOW+ITLINC
      ITUPR=ITUPR+ITUINC
      IF(ITLOW .LT. 0 .OR. ITLOW .GT. 255) GOTO 900 !OUTSIDE VALID RANGE
      IF(ITUPR .LT. 0 .OR. ITUPR .GT. 255) GOTO 900
      IF(ITLOW .GT. ITUPR) GOTO 900
C
C...   WRITE INFORMATION TO TERMINAL
C
      CALL SCROLL(1,24)
      CALL WRINF
      CALL CLRRLN(19,23)           !CLEAR LINES 19-23 (MENU LOCATION)
      CALL SCROLL(19,23)          !SET SCROLL REGION FOR 19-23
      GOTO 50
C
C...   RESET INITIAL THRESHOLD AND THEN RETURN
C
900    ITLOW=JTLW
      ITUPR=JTUPR
C
C...   PUT A <FF> IN FILE TO SEPERATE THIS RUN
C
      WRITE(2,910),12
910    FORMAT(' ',A1)

```

```

C...  WRITE SECOND PAGE OF INFO
      WRITE(2,915)
915   FORMAT('          AVERAGE
C          AVERAGE  AVERAGE'//
C          TOTAL    PERIMETER      A
CVERAGE    HORIZ    VERTICAL      AVERAGE
C          AVERAGE'//
C' THRESHOLD  PERIMETER  POINTS    TOTAL    PERI
CMETER      FERET      FERET      CUMULATIVE  IDD
C          ODmax'//
C' RANGE      POINTS    PER OBJECT  PERIMETER  PER O
C' OBJECT    DIAMETER  DIAMETER    IDD      PER OBJECT
C PER OBJECT'//)
      DO 930 J=1,NUMAVG
      WRITE(2,920), IAVG(J,1), IAVG(J,2), IAVG(J,3), GAVG(J,1), GAVG(J,7),
CGAVG(J,2), GAVG(J,3), GAVG(J,4), IAVG(J,4), GAVG(J,5), GAVG(J,6)
920   FORMAT(' ', I3, '- ', I3, 6X, I6, 5X, F10, 2, 4X, F10, 2, 6X, F7, 2, 6X, F7, 2,
C5X, F7, 2, 5X, I10, 3X, F12, 2, 4X, F7, 2)
930   CONTINUE
      WRITE(2,910), 12
      CLOSE(UNIT=2, DISPOSE='SAVE')      !CLOSE 'AVERAGE' FILE
C
999   CALL SCROLL(1,24)
      STATUS=1                          !RESET STATUS TO USER-INPUT
      CALL WRITF                          !WRITE INFO TO SCREEN
      CALL CLRLIN(18,23)
C
      RETURN
      END
C
C
C
C
C

```

```

SUBROUTINE CONTIG
  INCLUDE 'EDC_MSARCOM.FOR'
  BYTE OUTFIL(22),LOW(3),UPR(3),EXT(3)
  INTEGER*2 START(256),STOP(256),IADJPT(512),DROP(65536),MAP(65536)
  DATA EXT/'M','C','S'/

C
  UNIQUE=1
  IDROP=0

C
C... LOOP THROUGH ALL LINES
C
  DO 450 IROW=IRYTOP,IRYBOT
    TYPE 451,IROW
451  FORMAT(' ROW = 'I7)
  C
  C... SEARCH LINE FOR START AND STOP OF STRUCTURES ON LINE. THE ARRAYS
  C... START(*) AND STOP(*) CONTAIN THE STARTING AND STOPPING LOCATIONS
  C... OF ALL THE STRUCTURES ON THAT LINE. 'ISNUM' CONTAINS THE NUMBER OF
  C... STRUCTURES.
  C
    ISNUM=0
    JCOL=IRXTOP
90    IF (THRESH(JCOL,IROW) .EQ. 0) GOTO 140
    ISNUM=ISNUM+1
    START(ISNUM)=JCOL
100   JCOL=JCOL+1
    IF (JCOL .GT. IRXBOT) GOTO 110
    IF (THRESH(JCOL,IROW) .EQ. 1) GOTO 100
110   STOP(ISNUM)=JCOL-1
140   JCOL=JCOL+1
    IF (JCOL .LE. IRXBOT) GOTO 90

  C
  C... SEARCH EACH STRUCTURE FOR CONTIGUOUS STRUCTURES ABOVE
  C
    DO 400 ISTRCT=1,ISNUM
      IF (IROW .EQ. IRYTOP) GOTO 360          !TOP ROW

  C
  C... LOOP THROUGH STRUCTURE. SEARCH FOR ADJACENT STRUCTURES. IF FOUND,
  C... MAKE NOTE OF UNIQUE NUMBER OF ADJACENT STRUCTURE.
  C
    IADJ=0
    DO 300 JCOL=START(ISTRCT),STOP(ISTRCT)
      IF (THRESH(JCOL,IROW-1) .EQ. 0) GOTO 200
  C... ADJACENT POINT FOUND
      IF (IADJ .EQ. 0) GOTO 190
      DO 180 J=1,IADJ
        IF (THRESH(JCOL,IROW-1) .EQ. IADJPT(J)) GOTO 200 !UNIQUE # ALREADY FOUND
180     CONTINUE
190     IADJ=IADJ+1
        IADJPT(IADJ)=THRESH(JCOL,IROW-1)

  C
200   IF (IADJCT .EQ. 0) GOTO 300          !DONT CHECK FOR DIAGONALS
      IF (JCOL-1 .LT. IRXTOP) GOTO 230    !DONT CHECK OUTSIDE ROI
      IF (THRESH(JCOL-1,IROW-1) .EQ. 0) GOTO 230
  C... ADJACENT POINT FOUND
      IF (IADJ .EQ. 0) GOTO 220
      DO 210 J=1,IADJ
        IF (THRESH(JCOL-1,IROW-1) .EQ. IADJPT(J)) GOTO 230    !UNIQUE # ALREADY FOUND
210     CONTINUE
220     IADJ=IADJ+1

```

```

      IADJPT(IADJ)=THRESH(JCOL-1,IROW-1)
C
230  IF(JCOL+1.GT. IRXBOT) GOTO 300      !DONT CHECK OUTSIDE ROI
      IF(THRESH(JCOL+1,IROW-1).EQ. 0) GOTO 300
C...  ADJACENT POINT FOUND
      IF(IADJ.EQ. 0) GOTO 250
      DO 240 J=1,IADJ
      IF(THRESH(JCOL+1,IROW-1).EQ. IADJPT(J)) GOTO 300      !UNIQUE # ALREADY FOUND
240  CONTINUE
250  IADJ=IADJ+1
      IADJPT(IADJ)=THRESH(JCOL+1,IROW-1)
300  CONTINUE
C
      IF(IADJ.EQ. 0) GOTO 360      !NO ADJACENT POINTS. ASSIGN NEXT VALUE
      IF(IADJ.EQ. 1) GOTO 380      !ONE ADJACENT POINT
C
C...  MORE THAN ONE STRUCTURE IS ADJACENT. SET PRESENT STRUCTURE TO LOWEST
C...  ADJACENT STRUCTURE VALUE, THEN SEARCH FOR AND RESET ALL OTHER PREVIOUS
C...  STRUCTURE VALUES TO THAT SAME LOW VALUE. SEARCH BY GOING FROM PREVIOUS
C...  LINE TO THE FIRST LINE (BOTTOM TO TOP). WHEN STRUCTURE VALUE IS NO
C...  LONGER ENCOUNTERED, DROP IT FROM THE SEARCH
C
      IMIN=99999
      DO 310 J=1,IADJ      !SELECT MINIMUM STRUCTURE VALUE
      IF(IADJPT(J).LT. IMIN) IMIN=IADJPT(J)
310  CONTINUE
C
      DO 315 J=1,IADJ      !MAKE NOTE OF DROPPED VALUES
      IF(IADJPT(J).EQ. IMIN) GOTO 315
      IDROP=IDROP+1
      DROP(IDROP)=IADJPT(J)
315  CONTINUE
C
      DO 320 JCOL=START(ISTRCT),STOP(ISTRCT) !ASSIGN MIN VALUE TO 'THRESH'
      THRESH(JCOL,IROW)=IMIN
320  CONTINUE
C
      DO 350 K=1,IADJ
      IF(IADJPT(K).EQ. IMIN) GOTO 340      !DON'T DROP IMIN VALUE
      DO 340 I=IROW-1,IRYTOP,-1
      IFLG=0
      DO 330 J=IRXTOP,IRXBOT
      IF(THRESH(J,I).NE. IADJPT(K)) GOTO 330
      THRESH(J,I)=IMIN      !RESET VALUE
      IFLG=1
330  CONTINUE
      IF(IFLG.EQ. 0) GOTO 350      !DIDN'T FIND ON LINE. STOP SEARCH
340  CONTINUE
350  CONTINUE
      GOTO 400
C
C...  ASSIGN THE STRUCTURE NEXT UNIQUE VALUES. NO ADJACENT POINTS
C
360  DO 370 JCOL=START(ISTRCT),STOP(ISTRCT)
      THRESH(JCOL,IROW)=UNIQUE
370  CONTINUE
      UNIQUE=UNIQUE+1
      GOTO 400
C
C...  ASSIGN THE STRUCTURE THE VALUE OF THE SINGLE ADJACENT STRUCTURE

```

```

C
380      DO 390 JCOL=START(ISTRCT),STOP(ISTRCT)
          THRESH(JCOL, IROW)=IADJPT(1)
390      CONTINUE
          GOTO 400
C
400      CONTINUE
450      CONTINUE
C
C...    AT THIS POINT, THRESH CONTAINS UNIQUELY LABELLED MICROSTRUCTURES.
C...    HOWEVER, THE UNIQUE NUMBERS MUST BE COMPRESSED SINCE MANY VALUES
C...    WERE DROPPED. THE VALUES DROPPED ARE IN THE ARRAY 'DROP(IDROP)'.
C...    THIS IS ACCOMPLISHED BY CREATING A 'MAP', COMPRESSING THE MAP AND
C...    THEN USING THAT COMPRESSED MAP TO ALTER THE VALUES IN 'THRESH'
C
          UNIQUE=UNIQUE-1
C
          TYPE 452,UNIQUE
452      FORMAT(' NUMBER OF MSTs BEFORE COMPRESSION = 'I10//':')
          TYPE 454,UNIQUE, IDROP
454      FORMAT(' UNIQUE = 'I7' IDROP = 'I7)
          ICTMAP=0
C
          MAP(0)=0                                !CLEAR MAP
          DO 600 J=1,UNIQUE
600      MAP(J)=0
C
          DO 610 JD=1, IDROP                        !SETUP DROPPED POINTS
610      MAP(DROP(JD))=-1
C
          DO 620 J=1,UNIQUE                          !FILL IN MAP
          IF(MAP(J) .EQ. -1) THEN
              MAP(J)=0
              GOTO 620
          END IF
          ICTMAP=ICTMAP+1
          MAP(J)=ICTMAP
620      CONTINUE
          UNIQUE=ICTMAP
C
          DO 640 IROW=IRYTOP, IRYBOT                !APPLY MAP TO 'THRESH'
          DO 640 JCOL=IRXTOP, IRXBOT
          IF(THRESH(JCOL, IROW) .EQ. 0) GOTO 640
          THRESH(JCOL, IROW)=MAP(THRESH(JCOL, IROW))
640      CONTINUE
C
C...    WRITE TO OUTPUT FILE
C
          CALL CONFNM(OUTFIL,EXT)                  !GET OUTPUT FILENAME
          OPEN(UNIT=3,NAME=OUTFIL,TYPE='NEW')
          NCHANS=1                                !ONE CHANNEL OF DATA
          NHEADR=2                                !NUMBER OF HEADER LINES
          WRITE(3,650),NROWS,NCOLS,NCHANS,NHEADR
650      FORMAT(4I4)
C...    TYPE CHANNEL USED, ROI INFO, ADJACENCY FLAG AND THRESHOLD RANGE
          WRITE(3,670),ICHAN,IRXTOP,IRYTOP,IRXBOT,IRYBOT,IADJCT,ITLOW,ITUPR
670      FORMAT(B17)
C...    WRITE 'THRESH' ARRAY TO FILE
          DO 685 I=IRYTOP, IRYBOT
          WRITE(3,680), (THRESH(J, I), J=IRXTOP, IRXBOT)
685      CONTINUE

```

```
680  FORMAT(2216)
685  CONTINUE
      CLOSE(UNIT=3, DISPOSE='SAVE')
C
C
      RETURN
      END
C
C
C
C
```

```

SUBROUTINE STATIS(NUMAVG)
  INCLUDE 'EDC_MSARCOM.FOR'
  INTEGER*2 COORD(2,8), PERMCT, ORTHO, DIAG, ODmax, PERMC1, ICOUNT(8)
  INTEGER ODmaxS
  BYTE OUTFIL(22), EXT(3)
  DATA COORD/-1,-1,-1,0,-1,1,0,1,1,1,1,0,1,-1,0,-1/
  DATA EXT/'M','S','T'/
  SQRTWO=2**0.5

C...
C      OPEN FILE FOR MICROSTRUCTURE STATISTICS.  WRITE HEADER INFO
C
      CALL CONFMN(OUTFIL,EXT)          !GET OUTPUT FILENAME
      OPEN(UNIT=3,NAME=OUTFIL,TYPE='NEW')
      WRITE(3,10),OUTFIL,ITLOW,ITUPR
      FORMAT(' FILENAME = '22A1,T40' LOWER THRESHOLD = '13'  UPPER
10  C THRESHOLD = '13)
      WRITE(3,12),UNIQUE              !#MST,PIX SIZE, SIZE FLAG
12  FORMAT(' '15' UNIQUE MICROSTRUCTURES FOUND.')
      IF(IPIXFG.NE. 1) THEN
14          WRITE(3,14),SIZPIX,IPIXFG
          FORMAT(' ACTUAL PIXEL SIZE = 'F7.3' MICRONS.  FLAG = '12)
      ELSE
16          WRITE(3,16),SIZPIX,IPIXFG
          FORMAT(' DEFAULT PIXEL SIZE = 'F7.3' MICRONS.  FLAG = '12)
      END IF
      WRITE(3,22)
      WRITE(3,24)
      WRITE(3,26)
22  C      FORMAT('///'
      C
      C              MAXIMUM          WEIGHTED WEIGHTED
      C              UNIQUE          X          Y          X          Y
24  C      FORMAT(' PERIMETER PERIMETER HORIZONTAL VERTICAL
      C              AVERAGE')
26  C      FORMAT(' MICROSTRUCTURE CENTER CENTER CENTER CENTER
      C              AREA COUNT LENGTH RUN-LENGTH RUN-LENGTH IOD
      C              IOD ODmax')
C...
C      CLEAR SOME VALUES
C
      IAREAS=0
      XCEN1=0
      YCEN1=0
      XWCEN1=0
      YWCEN1=0
      XTCENT=0
      YTCENT=0
      PERMC1=0
      PERIM1=0
      IMAX_X_SUM=0
      IMAX_Y_SUM=0
      IODSUM=0
      ODmaxS=0

C...
C      CALCULATE STATISTICS FOR ALL MICROSTRUCTURES
C
      DO 900 IMST=1,UNIQUE

C...
C      SEARCH FOR FIRST OCCURANCE OF MICROSTRUCTURE #IMST
C
      DO 100 IROW=IRYTOP,IRYBOT

```

```

DO 100 ICOL=IRXTOP,IRXBOT
IF (THRESH(ICOL,IROW) .EQ. IMST) GOTO 110
CONTINUE
100
C
C... CALCULATE PERIMETER BE SEARCHING ADJACENT POINTS. IF IADJCT=0,
C... SEARCH #2,4,6,8. IF IADJCT=1, SEARCH #1,2,3,4,5,6,7,8.
C...
C...      1 2 3
C...      8 * 4
C...      7 6 5
C...
C... SINCE FIRST POINT WAS FOUND BY SEARCHING FROM THE LEFT (FROM #8)
C... WITH NOTHING IN THE PREVIOUS LINE, THE FIRST SEARCH POINT IS
C... #4. AFTER A POINT IS FOUND, SEARCH FOR THE NEXT POINT BY STARTING
C... AT (PREVIOUS SEARCH # + 4 + 1) OR (P.S.# - 4 + 1) IF P.S.# >= 4.
C
C
110 IPSRC=4
ISRC=IPSRC
PERMCT=0 !COUNT OF PERIMETER
ORTHO=0 !COUNT OF ORTHOGONAL POINTS
DIAG=0 !COUNT OF DIAGONAL POINTS
IYMIN=IROW !SET PARAMETERS FOR BLOB RANGE
IYMAX=IROW
IXMIN=ICOL
IXMAX=ICOL
IYST=IROW
IXST=ICOL
C
120 IF (IYST+COORD(1,ISRC) .LT. IRYTOP .OR. IYST+COORD(1,
CISRC) .GT. IRYBOT) GOTO 125 !CHECK IF OUT OF ROI
IF (IXST+COORD(2,ISRC) .LT. IRXTOP .OR. IXST+COORD(2,
CISRC) .GT. IRXBOT) GOTO 125
IF (THRESH (IXST+COORD(2,ISRC),IYST+COORD(1,ISRC)) .NE. 0) GOTO 160
125 ISRC=ISRC+1
IF (IADJCT .EQ. 0) ISRC=ISRC+1 !ONLY CHECK ORTHO POINTS
IF (ISRC .GT. 8) ISRC=ISRC-8
IF (ISRC .NE. IPSRC) GOTO 120
C... POINT IS ISOLATED. PERIMETER IS 1 POINT
PERMCT=1
GOTO 200
C
160 IF (JMOD (ISRC,2) .EQ. 0) THEN !DETERMINE IF ORTHO OR DIAO
ORTHO=ORTHO+1
ELSE
DIAG=DIAG+1
END IF
IYST=IYST+COORD(1,ISRC) !MOVE TO NEW POINT
IXST=IXST+COORD(2,ISRC)
IF (THRESH (IXST,IYST) .GT. 0) THEN
PERMCT=PERMCT+1 !COUNT 'NEW' POINT
THRESH (IXST,IYST) = -THRESH (IXST,IYST) !NEGATE POINT AS FLAG
END IF
C
C... AT THIS POINT, YOU MUST CHECK IF ALL SURROUNDING POINTS HAVE BEEN
C... COUNTED. SEARCH AND SEE IF ALL SURROUNDING POINT OF THE PRESENT
C... MICROSTRUCTURE HAVE BEEN SET TO -VALUE
C
IF (IYST .EQ. IROW .AND. IXST .EQ. ICOL) THEN

```

```

      IADJ_COUNT=0
      IADJ_NOT_COUNT=0
      ISKIP=2
      IF(IADJCT.EQ. 1) ISKIP=1
      DO 170 ICHK=ISKIP,8,ISKIP
      IF(IYST+COORD(1,ICLK) .LT. IRYTOP .OR. IYST+COORD(1,
C1CHK) .GT. IRYBOT) GOTO 170      !CHECK IF OUT OF ROI
      IF(IXST+COORD(2,ICLK) .LT. IRXTOP .OR. IXST+COORD(2,
C1CHK) .GT. IRXBOT) GOTO 170
      IF(THRESH(IXST+COORD(2,ICLK),IYST+COORD(1,ICLK)) .LT. 0) THEN
          IADJ_COUNT=IADJ_COUNT+1
          ICOUNT(IADJ_COUNT)=ICLK
      END IF
      IF(THRESH(IXST+COORD(2,ICLK),IYST+COORD(1,ICLK)) .GT. 0)
C IADJ_NOT_COUNT=IADJ_NOT_COUNT+1
170  CONTINUE
      IF(IADJ_NOT_COUNT.EQ. 0) GOTO 200      !BACK AT 1ST POINT
      IF(IADJ_COUNT.EQ. 1) GOTO 180      !CONTINUE SEARCH
      DO 172 JJ=1,IADJ_COUNT-1
      DO 172 KK=2,IADJ_COUNT
      IDIF=ABS(ICOUNT(JJ)-ICOUNT(KK))
      IF(IDIF.EQ. 1 .OR. IDIF.EQ. 7) GOTO 172
      GOTO 200      !SEPERATED BY TWO. BACK AT 1ST POINT
172  CONTINUE
      GOTO 180
      END IF

C
C
C      ISKIP=2
C      IF(IADJCT.EQ. 1) ISKIP=1
C      DO 170 ICHK=ISKIP,8,ISKIP
C      IF(IYST.LT. 60) GOTO 173
C      TYPE 172,IXST,IYST,ICLK,THRESH(IXST+COORD(2,ICLK),IYST+COORD(1,ICLK))
C172  FORMAT(' IXST = '15' IYST = '15' ICLK = '14' THRESH = '17)
C173  IF(THRESH(IXST+COORD(2,ICLK),IYST+COORD(1,ICLK)) .EQ. 0) THEN
C          IZRFLG=0
C          GOTO 170
C      END IF
C      IF(THRESH(IXST+COORD(2,ICLK),IYST+COORD(1,ICLK)) .LT. 0) IZRFLG=1
C      IF(THRESH(IXST+COORD(2,ICLK),IYST+COORD(1,ICLK)) .GT. 0
C      AND. IZRFLG.EQ. 0) GOTO 180
C170  CONTINUE
C      GOTO 200      !BACK AT 1ST POINT
C      END IF
C

C...  SET RANGE VALUES
180  IF(IYST.LT. IYMIN) IYMIN=IYST
      IF(IYST.GT. IYMAX) IYMAX=IYST
      IF(IXST.LT. IXMIN) IXMIN=IXST
      IF(IXST.GT. IXMAX) IXMAX=IXST
C...  CALCULATE WHERE TO LOOK NEXT
      IF(ISRC.LE. 4) THEN
          IPSRC=ISRC+4
      ELSE
          IPSRC=ISRC-4
      END IF
      IPSRC=IPSRC+2
      IF(IPSRC.GT. 8) IPSRC=IPSRC-8
      ISRC=IPSRC

```

```

GOTO 120
C
C... CALCULATE PERIMETER DISTANCE
C
200 PERIM=(ORTHO+(DIAG*SQRTWO))*SIZE1: !NORMALIZED TO MICRONS
C
C... SECOND PASS THROUGH RANGE THAT CONTAINS MICROSTRUCTURE.
C... DETERMINE MAX HORIZONTAL AND VERTICAL DIAMETERS
C
IMAX_Y=0
IMAX_X=0
DO 320 ICOL=IXMIN,IXMAX !CALCULATE Y-MAX DIAMETER
IST=0
DO 300 IROW=IYMIN,IYMAX
IF(IABS(THRESH(ICOL,IROW)) .EQ. IMST) THEN
IF(IST .EQ. 1) GOTO 290
IST=1
ISTLC=IROW
CONTINUE
290 ELSE
IF(IST .EQ. 0) GOTO 295
IF((IROW-1)-ISTLC .GT. IMAX_Y) IMAX_Y=IROW-1-ISTLC
IST=0
CONTINUE
295 END IF
300 CONTINUE
IF(IST .EQ. 0) GOTO 320
IF((IROW-1)-ISTLC .GT. IMAX_Y) IMAX_Y=IROW-1-ISTLC
CONTINUE
320
C
C... THIRD PASS THROUGH RANGE THAT CONTAINS MICROSTRUCTURE. DETERMINE
C... GEOMETRIC CENTER, AREA, IOD, ODmax, AND AVERAGE OPTICAL DENSITY.
C...
C... CLEAR ALL ARRAYS
IAREA=0
IYSUM=0
IXSUM=0
WYSUM=0.
WXSUM=0.
IOD=0
ODmax=-1
C
DO 420 IROW=IYMIN,IYMAX
IST=0
DO 400 ICOL=IXMIN,IXMAX
IF(IABS(THRESH(ICOL,IROW)) .EQ. IMST) THEN
IIMG=IMAGE(ICOL,IROW)
IF(IST .EQ. 1) GOTO 390
IST=1
ISTLC=ICOL
IAREA=IAREA+1
IYSUM=IYSUM+IROW
IXSUM=IXSUM+ICOL
IOD=IOD+IIMG
WYSUM=WYSUM+FLOAT(IROW*IIMG)
WXSUM=WXSUM+FLOAT(ICOL*IIMG)
IF(IIMG .GT. ODmax) ODmax=IIMG
390 ELSE
IF(IST .EQ. 0) GOTO 395
IF((ICOL-1)-ISTLC .GT. IMAX_X) IMAX_X=ICOL-1-ISTLC

```

```

395             1ST=0
               CONTINUE
END IF
400 CONTINUE
      IF(1ST.EQ. 0) GOTO 420
      IF((ICOL-1)-1STLC.GT. 1MAX_X) 1MAX_X=ICOL-1-1STLC
420 CONTINUE
      1MAX_Y=1MAX_Y+1
      1MAX_X=1MAX_X+1
      YCENT=FLOAT(1YSUM)/FLOAT(1AREA)           !GEOMETRIC CENTERS
      XCENT=FLOAT(1XSUM)/FLOAT(1AREA)
      YWCENT=WYSUM/FLOAT(1OD)
      XWCENT=WXSUM/FLOAT(1OD)
      AREA=FLOAT(1AREA)*(SIZPIX*SIZPIX)          !NORMALIZED AREA
      AVGO=FLOAT(1OD)/FLOAT(1AREA)               !AVERAGE OPTICAL DENSITY
      IF(1AREA.EQ. PERMCT) THEN                   !NO 'INTERIOR' POINTS
        PERMCT=-1
        PERIM=-1.0
      END IF
C
C... WRITE INFORMATION TO FILE
C
      WRITE(3,800), 1MST, XCENT, YCENT, XWCENT, YWCENT, AREA, PERMCT,
      CPERIM, 1MAX_X, 1MAX_Y, 1OD, AVGO, ODmax
800 FORMAT(' T6, 15, T18, F7. 2, T27, F7. 2, T36, F7. 2, T45, F7. 2, T53, F9. 2,
CT64, 15, T74, F8. 2, T89, 13, T102, 13, T109, 18, T118, F7. 2, T127, 13)
C
C... KEEP A RUNNING TOTAL OF ALL STATISTICS FOR THIS THRESHOLD RANGE.
C... THIS WILL BE WRITTEN TO AN OUTPUT FILE AT COMPLETION OF ALL PARTICLES
C
      XCEN1=XCEN1+XCENT           !X AND Y CENTERS
      YCEN1=YCEN1+YCENT
      XTCENT=XTCENT+(XCENT*1OD)
      YTCENT=YTCENT+(YCENT*1OD)
      XWCEN1=XWCEN1+XWCENT        !WEIGHTED X AND Y CENTERS
      YWCEN1=YWCEN1+YWCENT
      1AREAS=1AREAS+1AREA         !TOTAL AREA
      IF(PERMCT.LT. 0) GOTO 820
      PERMC1=PERMC1+PERMCT        !TOTAL PERIMETER COUNT
820 IF(PERIM.LT. 0) GOTO 830
      PERIM1=PERIM1+PERIM         !TOTAL PERIMETER LENGTH
830 1MAX_X_SUM=1MAX_X_SUM+1MAX_X  !TOTAL X,Y-RUN LENGTH
      1MAX_Y_SUM=1MAX_Y_SUM+1MAX_Y
      1ODSUM=1ODSUM+1OD           !TOTAL 1OD
      ODmaxS=ODmaxS+ODmax        !TOTAL ODmax
C
900 CONTINUE
      CLOSE(UNIT=3, DISPOSE='SAVE')
C
C... CALCULATE 'AVERAGE' OF STATISTICS AND WRITE TO UNIT 2
C
      IF(UNIQUE.LT. 1) THEN
        AUN0=0
        AOUT1=0
        AOUT2=0
        AOUT3=0
        AOUT6=0
        AOUT7=0
        GOTO 950
      END IF

```

```

AUNG=FLOAT(UNIQUE)
ADUT1=FLOAT(IAREAS)/AUNG
ADUT2=XCEN1/AUNG
ADUT3=YCEN1/AUNG
ADUT6=XTCENT/FLOAT(1ODSUM)
ADUT7=YTCENT/FLOAT(1ODSUM)
950  WRITE(2,920), ITLOW, ITUPR, IAREAS, UNIQUE, ADUT1, ADUT2, ADUT3,
      CADUT6, ADUT7
920  FORMAT(' ', 13, '-', 13, 6X, 16, 11X, 16, 11X, F9.2, 5X, F7.2, 6X, F7.2,
C6X, F7.2, 7X, F7.2)
C...  STORE STATS FOR SECOND SECTION OF OUTPUT
      NUMAVG=NUMAVG+1
      IAVG(NUMAVG, 1)=ITLOW
      IAVG(NUMAVG, 2)=ITUPR
      IF(UNIQUE .LT. 1) THEN
          IAVG(NUMAVG, 3)=0
          IAVG(NUMAVG, 4)=0
          GAVG(NUMAVG, 1)=0.
          GAVG(NUMAVG, 2)=0.
          GAVG(NUMAVG, 3)=0.
          GAVG(NUMAVG, 4)=0.
          GAVG(NUMAVG, 5)=0.
          GAVG(NUMAVG, 6)=0.
          GAVG(NUMAVG, 7)=0.
      ELSE
          IAVG(NUMAVG, 3)=PERMC1
          IAVG(NUMAVG, 4)=1ODSUM
          GAVG(NUMAVG, 1)=FLOAT(PERMC1)/AUNG
          GAVG(NUMAVG, 2)=PERIM1/AUNG
          GAVG(NUMAVG, 3)=FLOAT(IMAX_X_SUM)/AUNG
          GAVG(NUMAVG, 4)=FLOAT(IMAX_Y_SUM)/AUNG
          GAVG(NUMAVG, 5)=FLOAT(1ODSUM)/AUNG
          GAVG(NUMAVG, 6)=FLOAT(1ODmaxS)/AUNG
          GAVG(NUMAVG, 7)=PERIM1
      END IF
999  RETURN
      END
C
C
C
C
C

```

```

SUBROUTINE CONFNM(OUTFIL,EXT)
INCLUDE 'EDC_MSARCDM.FOR'
BYTE OUTFIL(22),LOW(3),UPR(3),EXT(3)

C      DO 700 JPOINT=1,14                                !COPY FILENAME
      IF(FNAME(JPOINT).EQ.' ') GOTO 710
      OUTFIL(JPOINT)=FNAME(JPOINT)
700    CONTINUE
710    ENCODE(3,3,LOW),ITLOW
      ENCODE(3,3,UPR),ITUPR
3      FORMAT(I3)
      OUTFIL(JPOINT)='_'
      OUTFIL(JPOINT+4)='_'
      DO 715 J=1,3
      IF(LOW(J).EQ.' ') LOW(J)='0'
      IF(UPR(J).EQ.' ') UPR(J)='0'
      OUTFIL(JPOINT+J)=LOW(J)
      OUTFIL(JPOINT+J+4)=UPR(J)
715    CONTINUE
      OUTFIL(JPOINT+8)='.'
      OUTFIL(JPOINT+9)=EXT(1)
      OUTFIL(JPOINT+10)=EXT(2)
      OUTFIL(JPOINT+11)=EXT(3)

C      RETURN
      END

C
C
C
C
C

```

```

SUBROUTINE WRTINF
  INCLUDE 'EDC_MSARCOM.FOR'
  BYTE BOLD(4), ATTOFF(4), REVERS(4), BLINK(4)
  DATA BOLD/27, 91, 49, 109/          !ESC [ 1 m
  DATA ATTOFF/27, 91, 48, 109/         !ESC [ 0 m
  DATA REVERS/27, 91, 55, 109/         !ESC [ 7 m
  DATA BLINK/27, 91, 53, 109/          !ESC [ 5 m

C
  ITOP=4
C
  CALL VTCURP(ITOP, 11)                  !FILENAME
  TYPE 100, FNAME
  CALL VTCURP(ITOP+4, 10)                 !LOWER THRESHOLD
  TYPE 140, ITLOW
  CALL VTCURP(ITOP+4, 27)                 !UPPER THRESHOLD
  TYPE 140, ITUPR
  CALL VTCURP(ITOP+6, 15)                 !LOWER_INCR
  TYPE 140, ITLIN
  CALL VTCURP(ITOP+6, 32)                 !UPPER_INCR
  TYPE 140, ITUINC
  CALL VTCURP(ITOP+4, 48)                 !ROWS
  TYPE 140, IRRDW
  CALL VTCURP(ITOP+4, 69)                 !COLUMNS
  TYPE 140, IRCOL
  CALL VTCURP(ITOP+5, 49)                 !X-TOP
  TYPE 140, IRXTOP
  CALL VTCURP(ITOP+5, 67)                 !Y-TOP
  TYPE 140, IRYTOP
  CALL VTCURP(ITOP+6, 52)                 !X-BOTTOM
  TYPE 140, IRXBOT
  CALL VTCURP(ITOP+6, 70)                 !Y-BOTTOM
  TYPE 140, IRYBOT
C
  CALL VTCURP(ITOP+5, 22)                 !THRESHOLD INCREMENT
  IF (ITINCR .EQ. 1) THEN
    TYPE 150
  ELSE
    TYPE 160
  END IF
C
  CALL VTCURP(ITOP+8, 21)                 !ADJACENCY CRITERIA
  IF (IADJCT .EQ. 1) THEN
    TYPE 170
  ELSE
    TYPE 180
  END IF
C
  CALL VTCURP(ITOP+9, 18)                 !PIXEL DIMENSION
  IF (IPIXFG .EQ. 1) THEN
    TYPE 190
  ELSE
    TYPE 200, SIZPIX
  END IF
C
  IF (STATUS .EQ. 1) THEN                 !PROGRAM STATUS
    TYPE 230, ATTOFF
    CALL VTCURP(ITOP+11, 17)
    TYPE 210
  ELSE
    TYPE 230, BOLD

```

```

                                TYPE 230, REVERS
                                TYPE 230, BLINK
                                CALL VTCURP(1TOP+11, 17)
                                TYPE 220
                                TYPE 230, ATTQFF
                                END IF
C
100  FORMAT('+', 14A1)
120  FORMAT('+', 11)
140  FORMAT('+', 13)
150  FORMAT('+ENABLED ')
160  FORMAT('+DISABLED ')
170  FORMAT('+ORTHOGONAL/DIAGONAL ')
180  FORMAT('+ORTHOGONAL ')
190  FORMAT('+++++')
200  FORMAT('+', F7.3)
210  FORMAT('+USER-INPUT ')
220  FORMAT('+EXECUTING ')
230  FORMAT('+', 4A1)
C
C
                                RETURN
                                END
C
C
C
C

```

```

SUBROUTINE VTCURP(LINE,COL)
C
C...  PLACE CURSOR AT LINE & COL OF SCREEN
C
      IMPLICIT INTEGER*2(A-Z)
      BYTE CP(8)    !HORIZONTAL & VERTICAL CURSOR POSITION
      DATA CP/27,91,48,48,59,48,48,72/ !ESC [ 0 0 ; 0 0 H
      ENCODE(2,2,CP(3))LINE
      ENCODE(2,2,CP(6))COL
      FORMAT(12)
2     IF (CP(3) .EQ. ' ') CP(3) = '0'
      IF (CP(6) .EQ. ' ') CP(6) = '0'
      TYPE 10,CP
10    FORMAT('*,BA1)
C
      RETURN
      END
C
C
C
C
C

```

```

SUBROUTINE CLRTXT
  BYTE COMAND(4)
  DATA COMAND/27,91,50,74/

C
100  TYPE 100,COMAND
      FORMAT('+'4A1)
      RETURN
      END

C
C
C
C
C
      SUBROUTINE CLRLIN(ITOP,IBOT)
      BYTE MOVE(7),ERASE(4)
      DATA MOVE/27,91,0,0,59,49,72/          !ESC [ P1 , 1 H
      DATA ERASE/27,91,50,75/                !ESC [ 2 K

C
      DO 200 J=ITOP,IBOT
      ILINE=J
      MOVE(3)='0'
      IF(ILINE.GT. 9) MOVE(3)=(ILINE/10)+48
      MOVE(4)=JMOD(ILINE,10)+48
      TYPE 100,MOVE          !MOVE CURSOR TO LINE
      TYPE 100,ERASE         !ERASE ALL OF LINE
100  FORMAT('+'7A1)
200  CONTINUE
      RETURN
      END

C
C
C
C
C
      SUBROUTINE CLRSCR
      BYTE INIT(3),END(2)
      DATA INIT/27,80,112/          !ESC P p
      DATA END/27,92/              !ESC \

C
      TYPE 100,INIT,END
100  FORMAT('+'3A1,'S(E)',2A1)
      RETURN
      END

C
C
C
C
C

```

```

C      SUBROUTINE VTDBLE(TEXT,NCHARS,LINE,COL,VID)
C...
C      DOUBLE WIDTH VT TERMINAL TEXT PROGRAM
C
C      CALLING VARIABLES:
C
C      TEXT      - input string of text
C      NCHARS    - input number of characters of the text [0-40]
C                  if nchars = 0, simply set the existing text to double width
C      LINE      - input line of text [ 1 - 23 ]
C      COL       - input starting column of text [ 1 - 40 ]
C      VID       - normal (0) or reverse (1) video
C
C      IMPLICIT INTEGER*2(A-Z)
C      BYTE DW(3),TEXT(40),CP(8),RV(4),NV(3)
C      DATA CP/27,91,48,48,59,48,48,72/ !CURSOR POSITION
C      DATA DW/27,35,54/ !DOUBLE WIDTH SINGLE HEIGHT
C      DATA RV/27,91,55,109/ !ESC [ 7 m REVERSE VIDEO
C      DATA NV/27,91,109/ !ESC [ m NORMAL VIDEO
C      IF (LINE .LT. 1) LINE = 1 !VT SCREEN LINES FROM 1
C      IF (LINE .GT. 23) LINE = 23 ! TO 23
C      IF (COL .LT. 1) COL = 1
C      IF (COL .GT. 40) COL = 40
C      IF (NCHARS + COL .GT. 41) NCHARS = 41 - COL
C      ENCODE(2,2,CP(3)) LINE
C      ENCODE(2,2,CP(6)) COL
C      2 FORMAT(I2)
C      IF (CP(3) .EQ. ' ') CP(3) = '0' !ZERO FILL 1 DIGIT NUMBERS
C      IF (CP(6) .EQ. ' ') CP(6) = '0'
C      IF (NCHARS .LE. 0) THEN !make existing text double width
C          TYPE 10,CP,DW
C          RETURN
C      end if
C      if (vid .eq. 1) then
C          type 10,Cp,Dw,Rv,(Text(i),i=1,Nchars),Nv !reverse video text
C      else
C          type 10,CP,DW,NV,(TEXT(I),I=1,NCHARS) !normal video text
C      end if
C      10 FORMAT('s',58A1)
C      RETURN
C      END
C
C
C
C

```

```
C
C
C ASYNCHRONOUS KEYBOARD INPUT USING SYSTEM SERVICES QIO
C WITH WAIT FOR EVENT FLAG TO PREVENT POLLING.
C .....
C
C FUNCTION IESC ( CHAR )
C
C IMPLICIT INTEGER*2 ( A - Z )
C LOGICAL*1 FLAG
C INTEGER*2 IESC_EFN
C INTEGER*2 IOSB(4)
C INTEGER*2 IESC_CHAN
C INTEGER*2 NEXT(4)
C BYTE FIRST_CHAR
C BYTE LAST_TWO_CHARS(2)
C
C CHARACTER*2 TERMINAL / 'TT' /
C
C DATA FLAG / .FALSE. /
C
C INCLUDE '( $IODEF )' !DEFINE SYMBOLIC NAMES FOR I/O FUNCTION
C CODE IO$READVBLK AND IO$M_NOECHO, THE
C CODE MODIFIER.
C .....
C
C CHAR = 0
C IESC = 0
C IESC_EFN = 2
C IF (FLAG) GO TO 10
C ISTAT = SYS$ASSIGN (TERMINAL,IESC_CHAN,,)
C IF (.NOT. ISTAT ) CALL LIB$STOP ( %VAL (ISTAT ) )
C FLAG = .TRUE.
C IO_FUNC = IO$_READVBLK .OR. IO$_M_NOECHO .OR. IO$_M_ESCAPE
C IO_TEST = IO$_SENSEMODE .OR. IO$_M_TYPEAHD CNT
C ISTAT = SYS$QIOW(%VAL(IESC_EFN), , ! READ A CHARACTER
C 1 %VAL(IESC_CHAN), ! WHEN KEYS.
C 2 %VAL(IO_FUNC ), !
C 3 %REF(IOSB ), !
C 4 , !
C 5 , !
C 6 %REF(FIRST_CHAR), !
C 7 %VAL(1), !
C 8 , !
C 9 , !
C 1 , ) !
C
C IF (.NOT. ISTAT ) CALL LIB$STOP (%VAL(ISTAT)) ! STOP IF QIO NOT SUCCESS
C CHAR = FIRST_CHAR
C ISTAT = SYS$QIOW(%VAL(IESC_EFN), !HOW MANY MORE CHARACTERS
C 1 %VAL(IESC_CHAN),
C 2 %VAL(IO_TEST),
C 3 %REF( IOSB ),
```

```

4      .
5      .
6      . %REF(NEXT),.... )
C
C      IF( NEXT(1) .EQ. 0 ) RETURN
C
C      ISTAT = SYS$QIOH(%VAL(IESC_EFN),          !GET THE REST
1      . %VAL(IESC_CHAN),
2      . %VAL(IESC_FUNC),
3      . %REF(      IOSB),
4      .
5      . %REF(LAST_TWO_CHARS),%VAL(2),.... )
6
C
C      CHAR = LAST_TWO_CHARS(2)
C      IESC = 1
C
C      WRITE (6,100) IESC,IDUM
CC100  FORMAT (6X,'BYTE COUNT: ' ,I4, ' FIRST CHARACTER: ',A1 )
C
C      20  RETURN
C      END
C
C
C
C
C

```

```

SUBROUTINE SCROLL(ITOP,IBOT)
C
C... SET SCROLL REGION TO ONE DEFINED BY ITOP,IBOT. THEN MOVE
C... CURSOR TO 'HOME', TOP LINE OF SCROLL REGION.
C
      BYTE INIT(2),STRING(6)
      DATA INIT/27,91/           !ESC [
C
      STRING(1)=(ITOP/10)+48
      STRING(2)=(ITOP-((STRING(1)-48)*10))+48
      STRING(3)=' '
      STRING(4)=(IBOT/10)+48
      STRING(5)=(IBOT-((STRING(4)-48)*10))+48
      STRING(6)='r'
C
      TYPE 10,"33,(STRING(J),J=1,6)
10    FORMAT('+',A1,'L',6A1)
      TYPE 20,INIT
20    FORMAT('+',2A1,'?6h')
      RETURN
      END

```